

IDENTIFYN RESOLUTION STEPDOWN

PARP1

MICROSCOPY EVIDENCE ASSESSMENT

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532

WESTERN BLOT -> WIDEFIELD -> WIDEFIELD SIM / APOTOME -> CONFOCAL -> AIRYSCAN
SIM -> SIM²

7


ORDERED EVIDENCE TIERS

0

SUPPORTING CONTROL


PENDING

SCIENTIST REVIEW



Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn™ SRM Alexa Fluor™ 532, IgG - BLR018M
Cat ID# - DC 000053

Panel A shows a three-dimensional volume projection of PARP1 detected with the Alexa Fluor™ 532 direct conjugate (white) and DAPI (blue). A white box denotes the region of interest (ROI) analyzed in Panel B, which presents a three-dimensional surface projection of the same PARP1 assemblies. At SIM² lateral resolution (~40–60 nm), PARP1 displays nuclear localization and aggregate morphologies comparable to those observed by lattice SIM, consistent with the relatively large size of PARP1 assemblies. Structural differentiation between SIM and SIM² is therefore limited at this stage; further resolution of PARP1 sub-units and higher-order organization requires single-molecule localization microscopy.



REPRESENTATIVE HIGHEST-AVAILABLE VIEW / SIM²

DC-000053 | PRODUCT-LEVEL EVIDENCE DOSSIER

SOURCE-GROUNDED MICROSCOPY + BIOLOGICAL REVIEW / IDENTIFYN EVIDENCE PRESERVED

CRO CAPABILITY DEMONSTRATED

Identifyn performed the experimental and imaging work documented below. BENNETT converts that operating history into a governed, source-bounded scientific record.

- Source-documented sample preparation and antibody/reagent workflow
- Western blot imaging
- Widefield fluorescence microscopy
- Widefield structured illumination / Apotome optical sectioning
- Confocal microscopy
- Airyscan acquisition and reconstruction workflow
- Structured illumination microscopy
- SIM² reconstruction workflow
- Treatment-response experimental workflow
- Z-stack / 3D acquisition
- Multi-channel imaging

SCIENTIFIC QUESTION

What source-bounded experimental and microscopy evidence is preserved for PARP1, and what is the strongest interpretation allowed by the available modalities and controls?

EVIDENCE AVAILABLE

Western blot -> Widefield -> Widefield SIM — Apotome -> Confocal -> Airyscan -> SIM -> SIM²

BENNETT ASSESSMENT

The preserved DC-000053 record contains the applicable NIR sequence: Western blot; Widefield; Widefield SIM — Apotome; Confocal; Airyscan; SIM; SIM². The historical boundary reflects the validated application available for the specific conjugate, deployed platform and workflow, not a universal scientific impossibility.

BENNETT / MICROSCOPY AND BIOLOGICAL EVIDENCE REVIEW

SOURCE-GROUNDED RETROSPECTIVE REVIEW / PENDING SCIENTIST REVIEW

SOURCE STATE	SOURCE HASH VERIFIED / EVENT RECORDED
BENNETT EVIDENCE REVIEW	COMPLETED / SOURCE-BOUNDED
NATIVE IMAGE REPROCESSING	NOT PERFORMED
LITERATURE SOURCES	INDEPENDENTLY VERIFIED
SCIENTIST REVIEW	PENDING
PUBLICATION	NOT AUTHORIZED

OVERALL CALL / NIR-BOUNDED PARTIAL STEPDOWN

The preserved DC-000053 record contains the applicable NIR sequence: Western blot; Widefield; Widefield SIM — Apotome; Confocal; Airyscan; SIM; SIM². The historical boundary reflects the validated application available for the specific conjugate, deployed platform and workflow, not a universal scientific impossibility.

SETTINGS ASSESSMENT

The record preserves 7 ordered evidence tier(s). Structured source metadata values are reported by stage and remain tied to the original method objects.

BIOLOGICAL SUPPORT

The record supports retrospective, source-bounded microscopy review. Target identity, specificity, treatment effect, binding, interaction and mechanism are not inferred from presentation images alone.

CONTROL ASSESSMENT

No separate control record is preserved; none is inferred from method or template language.

CLAIM CEILING

This draft supports historical capability documentation, source registration, modality ordering, settings review and bounded interpretation. It does not independently validate antibody specificity, quantify a population response, establish binding, interaction, causation or mechanism, or authorize publication.

PROCESSING DISCLOSURE

Retrospective BENNETT architectural evaluation of preserved Identifyn source evidence. BENNETT did not perform native-file reprocessing, new deterministic measurement or the historical image post-processing represented in this record.

REAGENTS AND ACQUIRED CHANNELS

Preparation reagents are reported separately from channels actually documented in the acquisition settings. Fluorophores mentioned in staining are not automatically treated as acquired channels.

EVIDENCE OBJECT	REAGENTS / FLUOROPHORES PRESENT	CHANNELS ACTUALLY ACQUIRED
Western blot	532	Not deterministically stated in the acquisition block
Widefield	405/DAPI; 532; 555; 647	3; AF532-49014; 50 Cy 5; 49 DAPI; 515 - 545; 625 - 655; 335 - 383; 555 - 595
Widefield SIM — Apotome	405/DAPI; 532; 555; 647	3; AF532-49014; 50 Cy 5; 49 DAPI; 515 - 545; 625 - 655; 335 - 383; 555 - 595
Confocal	405/DAPI; 488; 532; 647	3; None; 488nm @0.8%; 639nm @0.10%; 405nm @0.2%; 528; 653; 353
Airyscan	405/DAPI; 488; 532; 647; 750	3; None; 488nm @1.20%; 639nm @0.10%; 405nm @0.2%; 528; 653; 353
SIM	405/DAPI; 488; 532; 647	3; SR HE LBF 405/488/561/642; 200 - 405, 200 - 488, 200 - 561, 200 - 642; 419 - 470, 503 - 546, 576 - 617, 657 - 800; T2-Axiocam 820 #2-SR; T1-Axiocam 820 -SR; T3-Axiocam 820 #2-SR; 488
SIM ²	405/DAPI; 488; 532; 647	3; SR HE LBF 405/488/561/642; 200 - 405, 200 - 488, 200 - 561, 200 - 642; 419 - 470, 503 - 546, 576 - 617, 657 - 800; T2-Axiocam 820 #2-SR; T1-Axiocam 820 -SR; T3-Axiocam 820 #2-SR; 488

MODALITY ASSESSMENT / PART 1 OF 7

Western blot / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: Azure Sapphire FL. Objective: not explicitly stated. Platform: Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s). Structured source metadata values: 11.

VISIBLE OBSERVATION

A preserved presentation image is available for Western blot. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Western blot, documents platform context as Azure Sapphire FL, and records detected dye/channel descriptors as 532.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

Band position and displayed intensity do not independently establish analytical specificity, target identity, linearity, or treatment effect.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Western blot evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 2 OF 7

Widefield / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Axio Imager.M1. Objective: Plan-Apochromat 63x/1.40 Oil M27. Platform: Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705 mono, and X-Cite Xylis Lamp, was used with the following parameters for each dye: Acquisition: Channels: 3; Scaling (Per Pixel): 0.110µm X 0.110µm; Image size (Pixels): 725 X 644; Image Size (Pixels): 725 X 644; Bit Depth: 14 Bit. Detector/camera: Imaging Device: Axiocam 705. Timing: Exposure Time: 250 ms; 25 ms. Illumination: Light source: X-Cite Xylis Lamp Intensity @25.00%; X-Cite Xylis Lamp Intensity @10.00%. Optical/scan: Effective NA: 1.4. Structured source metadata values: 37.

VISIBLE OBSERVATION

A preserved presentation image is available for Widefield. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Widefield, documents platform context as ZEISS Axio Imager.M1, and records detected dye/channel descriptors as 405/DAPI; 532; 555; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

A presentation image supports gross localization and source-bounded co-occurrence, not direct binding, molecular colocalization, or native quantitative measurement.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Widefield evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 3 OF 7

Widefield SIM — Apotome / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Axio Observer.Z1/7 Apotome. Objective: EC Plan-Neofluar 63X/1.25 oil M27. Platform: Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:. Acquisition: Z-Stack: 98 Slices (9.7 μm); Channels: 3; Scaling (Per Pixel): 0.143 μm X 0.143 μm X 0.100 μm ; Image size (Pixels): 535 X 417; Image Size (Pixels): 535 X 417; Bit Depth: 14 Bit. Detector/camera: Imaging Device: Axiocam 807 Mono. Timing: Exposure Time: 200 ms; 25 ms. Illumination: Light source: X-Cite Xylis Lamp Intensity @20.00%; X-Cite Xylis Lamp Intensity @10.00%. Optical/scan: Effective NA: 1.4; 1.25. Processing: Projection: View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 56.

VISIBLE OBSERVATION

A preserved presentation image is available for Widefield SIM — Apotome. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Widefield SIM — Apotome, documents platform context as ZEISS Axio Observer.Z1/7 Apotome, and records detected dye/channel descriptors as 405/DAPI; 532; 555; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

Improved optical sectioning does not by itself establish reconstructed super-resolution or molecular interaction.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Widefield SIM — Apotome evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 4 OF 7

Confocal / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS LSM 980. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:. Acquisition: Z-Stack: 121 Slices (3.6 μm); Channels: 3; Scaling (Per Pixel): 0.059 μm X 0.059 μm X 0.030 μm ; Image size (Pixels): 1121 X 1121; Image Size (Pixels): 1121 X 1121; Bit Depth: 16 Bit. Detector/camera: Imaging Device: GaAsP-Pmt3; GaAsP-Pmt3; Detector Type: GaAsP-PMT; GaAsP-PMT. Timing: Pixel Time: 0.30 μs ; Frame Time: 3.69s. Illumination: Laser Wavelength: 488nm @0.8%; 639nm @0.10%. Optical/scan: Pinhole: 1.00AU / 53 μm ; 1.00AU / 66 μm ; Scan Zoom: 2.0; LSM Scan Speed: 12; Scan Direction: Bidirectional; Averaging: 8; Effective NA: 1.4. Processing: Projection: View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black; LSM Plus Mode: Filter: 6.8 (3D, Auto); Filter: 7.5 (3D, Auto). Structured source metadata values: 59.

VISIBLE OBSERVATION

A preserved presentation image is available for Confocal. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Confocal, documents platform context as ZEISS LSM 980, and records detected dye/channel descriptors as 405/DAPI; 488; 532; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

Optical sectioning and source-bounded 3D presentation do not establish binding, mechanism, or population-level effect.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Confocal evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 5 OF 7

Airyscan / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS LSM 980. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:. Acquisition: Z-Stack: 141 Slices (4.2 μm); Channels: 3; Scaling (Per Pixel): 35.30 nm X 35.30 nm X 30.00 nm; Image size (Pixels): 607 X 798; 241 X 178; Image Size (Pixels): 607 X 798; 241 X 178; Bit Depth: 16 Bit. Detector/camera: Imaging Device: Airyscan; Detector Type: GaAsP-PMT; GaAsP-PMT. Timing: Pixel Time: 0.66 μs ; Frame Time: 7.82s. Illumination: Laser Wavelength: 488nm @1.20%; 639nm @0.10%. Optical/scan: Pinhole: 5.00AU / 265 μm ; 5.00AU / 328 μm ; Scan Zoom: 2.4; LSM Scan Speed: 7; Scan Direction: Bidirectional; Averaging: 4; Effective NA: 1.4. Processing: Projection: View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 59.

VISIBLE OBSERVATION

A preserved presentation image is available for Airyscan. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Airyscan, documents platform context as ZEISS LSM 980, and records detected dye/channel descriptors as 405/DAPI; 488; 532; 647; 750.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Displayed feature separation remains reconstruction dependent; no unsupported numerical resolution or molecular architecture claim is permitted.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Airyscan evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 6 OF 7

SIM / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Elyra. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye: Acquisition: Z-Stack: 142 Slices (4.23 μ m); Channels: 3; Scaling (Per Pixel): 21.11 nm X 21.11 nm X 30.00 nm; Image size (Pixels): 416 X 462; 184 X 164; Image Size (Pixels): 416 X 462; 184 X 164; Bit Depth: 16 Bit. Detector/camera: Imaging Device: AxioCam 820 #2; AxioCam 820. Timing: Exposure Time: 150 ms; 100 ms; Frame Time: 1.59 s; 1.33 s. Illumination: Laser Wavelength: 488 nm @4.0%; 640 nm @2.0%; Illumination Wavelength: 488; 640. Optical/scan: Scan Direction: Unidirectional; Effective NA: 1.4. Processing: Algorithm: SIM; Process Sampling: N/A; Result Sampling: 4; Projection: View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); View Angle 35°, Scale Z 35%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 65.

VISIBLE OBSERVATION

A preserved presentation image is available for SIM. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as SIM, documents platform context as ZEISS Elyra, and records detected dye/channel descriptors as 405/DAPI; 488; 532; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Reconstruction-dependent structural separation requires native reconstruction QC before numerical resolution or molecular dimensions can be claimed.

CLAIM CEILING

This tier supports a source-bounded statement that the historical SIM evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 7 OF 7

SIM² / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Elyra. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye: Acquisition: Z-Stack: 110 Slices (3.27 μ m); Channels: 3; Scaling (Per Pixel): 21.11 nm X 21.11 nm X 30.00 nm; Image size (Pixels): 441 X 358; 192 X 218; Image Size (Pixels): 441 X 358; 192 X 218; Bit Depth: 16 Bit. Detector/camera: Imaging Device: Axiocam 820 #2; Axiocam 820. Timing: Exposure Time: 150 ms; 100 ms; Frame Time: 1.59 s; 1.33 s. Illumination: Laser Wavelength: 488 nm @4.0%; 640 nm @2.0%; Illumination Wavelength: 488; 640. Optical/scan: Scan Direction: Unidirectional; Effective NA: 1.4. Processing: Algorithm: SIM2; Process Sampling: 4; Result Sampling: 4; Projection: View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 67.

VISIBLE OBSERVATION

A preserved presentation image is available for SIM². BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as SIM², documents platform context as ZEISS Elyra, and records detected dye/channel descriptors as 405/DAPI; 488; 532; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Further reconstruction-dependent segmentation does not establish molecular dimensions, structure identity, or mechanism without independent support.

CLAIM CEILING

This tier supports a source-bounded statement that the historical SIM² evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

CROSS-MODALITY SYNTHESIS

The preserved DC-000053 record contains the applicable NIR sequence: Western blot; Widefield; Widefield SIM — Apotome; Confocal; Airyscan; SIM; SIM². The historical boundary reflects the validated application available for the specific conjugate, deployed platform and workflow, not a universal scientific impossibility.

Different acquisitions are not treated as a quantitative intensity ladder. Any continuity statement remains qualitative and source bounded.

BENNETT LITERATURE CORRELATION — SOURCES INDEPENDENTLY VERIFIED

BENNETT uses independently verified primary scientific sources to test whether the interpretation is scientifically consistent and properly bounded. Literature correlation does not retroactively validate the Identifyn dataset or replace scientist review.

- Antibody and microscopy evidence require orthogonal controls and modality-appropriate interpretation; presentation images alone do not establish analytical specificity.
- Advanced-resolution presentation depends on acquisition, reconstruction or localization quality control before numerical resolution, molecular dimensions or cluster identity can be claimed.

REFERENCES

1. Uhlen M, et al. A proposal for validation of antibodies. Nat Methods. 2016;13:823-827. doi:10.1038/nmeth.3995.
2. Gustafsson MGL. Surpassing the lateral resolution limit by a factor of two using structured illumination microscopy. J Microsc. 2000;198:82-87. doi:10.1046/j.1365-2818.2000.00710.x.
3. Sage D, et al. Quantitative evaluation of software packages for single-molecule localization microscopy. Nat Methods. 2015;12:717-724. doi:10.1038/nmeth.3442.

RELEASE BOUNDARY

This assessment remains a draft pending scientist review. No publication is authorized. Any future native-file analysis must enter BENNETT as a new evidence snapshot with routed engine authority, native-data QA, methods, limitations, and an independently reviewed claim ceiling.

SOURCE RECONCILIATION / QA / PROVENANCE

SOURCE STATUS - PRESERVED / AUTO-RECONCILED

Historical Identifyn source material was subjected to BENNETT retrospective QA. Original source state and correction history remain preserved in provenance.

HIGH CONFIDENCE - AUTO-RECONCILED

Airyscan / Scaling (Per Pixel): 35.30 nm X 35.30 nm X 30.00 nm -> 0.03530 μ m X 0.03530 μ m X 0.03000 μ m. Both numerical magnitude and unit were normalized, then checked against the reported field dimensions; no unit label was changed alone. Supporting evidence: Image Size (Pixels)=607 X 798; Image Size (Scaled)=21.42 μ m X 28.14 μ m; PASS - INTERNALLY CONSISTENT. Documented sampling reproduces the reported physical field dimensions. Maximum field-dimension error 0.10%.

HIGH CONFIDENCE - AUTO-RECONCILED

SIM / Scaling (Per Pixel): 21.11 nm X 21.11 nm X 30.00 nm -> 0.02111 μ m X 0.02111 μ m X 0.03000 μ m. Both numerical magnitude and unit were normalized, then checked against the reported field dimensions; no unit label was changed alone. Supporting evidence: Image Size (Pixels)=416 X 462; Image Size (Scaled)=8.78 μ m X 9.75 μ m; PASS - INTERNALLY CONSISTENT. Documented sampling reproduces the reported physical field dimensions. Maximum field-dimension error 0.03%.

HIGH CONFIDENCE - AUTO-RECONCILED

SIM² / Scaling (Per Pixel): 21.11 nm X 21.11 nm X 30.00 nm -> 0.02111 μ m X 0.02111 μ m X 0.03000 μ m. Both numerical magnitude and unit were normalized, then checked against the reported field dimensions; no unit label was changed alone. Supporting evidence: Image Size (Pixels)=441 X 358; Image Size (Scaled)=9.31 μ m X 7.56 μ m; PASS - INTERNALLY CONSISTENT. Documented sampling reproduces the reported physical field dimensions. Maximum field-dimension error 0.03%.

HIGH CONFIDENCE - SCIENTIST-AUTHORIZED CORRECTION

Cross-tier / treatment concentration: 200 μ M; 2 μ M; 2 μ M; 647 nm -> Etoposide = 2 μ M. The governing scientist identified every 200 μ M occurrence as a technician transcription error. The canonical historical concentration is 2 μ M. Supporting evidence: Scientist-authorized correction supplied for the Identifyn historical archive; original technician text remains preserved unchanged. Authority: Dr. Brian T. Bennett. Preservation: Original source text and hashes remain unchanged; only the derivative canonical field is corrected.

UNRESOLVED - PRESERVED AS CONFLICT

Cross-tier / treatment duration: {"Western blot": ["20 hours"], "Widefield": ["24 hours"], "Widefield SIM — Apotome": ["24 hours"], "Confocal": ["24 hours"], "Airyscan": ["24 hours"], "SIM": ["24 hours"], "SIM²": ["24 hours"]} -> Preserved as modality-specific historical duration values. Genuinely different treatment durations are present across evidence stages, and no source or scientist correction resolves them. Supporting evidence: Only duration values from explicit etoposide-treatment sentences were classified; concentration, wavelength, resolution and pixel-size values were excluded.

BENNETT EVIDENCE REVIEW COMPLETE

INTERPRETATION TO EVIDENCE

BENNETT has completed the bounded retrospective assessment.
The following pages switch ownership to the preserved Identifyn record.



IDENTIFYN PRESERVED SOURCE EVIDENCE

ORIGINAL RECORD BEGINS

Images, source methods, acquisition settings and provenance follow in the locked
Stepdown order. Scientific wording is preserved; missing modalities are not invented.
No separate control record is present in the preserved product record.

PRODUCT-LEVEL STEPDOWN MAP


This dossier preserves the complete available evidence progression for DC-000053. Each tier remains tied to its original image, source method, instrument context and cryptographic identity. No separate control record is present in the preserved product record, and none is inferred.

CANONICAL STEPDOWN	EVIDENCE STAGE	INSTRUMENT	STATUS
1	Western blot	Azure Sapphire FL	PRESERVED
2	Widefield	ZEISS Axio Imager.M1	PRESERVED
3	Widefield SIM - Apotome	ZEISS Axio Observer.Z1/7 Apotome	PRESERVED
4	Confocal	ZEISS LSM 980	PRESERVED
5	Airyscan	ZEISS LSM 980	PRESERVED
6	SIM	ZEISS Elyra	PRESERVED
7	SIM ²	ZEISS Elyra	PRESERVED

INTERPRETIVE BOUNDARY

Ordering evidence by resolution tier does not make the modalities interchangeable. Each stage retains its acquisition environment and scientific limitations. This pilot organizes the record; it does not synthesize new biological claims.

ORIGINAL IMAGE EVIDENCE / WESTERN BLOT



EVIDENCE TIER	Western blot
INSTRUMENT	Azure Sapphire FL
CELL MODEL	HeLa
DETECTED DYES	532
METADATA VALUES	11

STRUCTURED ACQUISITION METADATA / WESTERN BLOT

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)
Scanner	Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)
Detection mode	Fluorescence
Laser	532nm
Emission Filter	624±40nm
Detector	PMT
Intensity	2
Pixel size	100 µm
Scan Speed	High
Quality	High
Focus	Membrane (0 µm in Z-axis)
Image	K. Small

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Expected MW: 116kD

HeLa cells were treated with 200 µM etoposide for 20 hours to induce DNA damage. Cells were then homogenized using an ultrasonication probe and RIPA lysis buffer, and then centrifuged to pellet the cell debris. The lysate supernatant (100 µg) was mixed with NuPAGE™ LDS Sample Buffer (4X), heated at 80 °C for 10 min, and subjected to electrophoresis using NuPAGE™ 4 to 12% Bis-Tris, 1.0-1.5 mm, Mini Protein Gels in a Mini Gel Tank. The Thermo Fisher's PageRuler™ Prestained Protein Ladder, 10 to 180 kDa, was used as a molecular weight marker.

Post electrophoresis, the proteins were transferred to a polyvinylidene difluoride (PVDF) membrane using the iBlot™ 3 Starter Kit PVDF. The transfer membrane was blocked with Blocker™ Casein for 1 hour and then washed five times with PBS. The blocked membrane was incubated with Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate for 2 hours, followed by five washes in 1X PBS. The membrane was dried for 10 minutes at room temperature before imaging.

Scanner = Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)

Scan Settings

Detection mode = Fluorescence
Laser = 532nm
Emission Filter = 624±40nm
Detector = PMT
Intensity = 2
Pixel size = 100 µm
Scan Speed = High
Quality = High
Focus = Membrane (0 µm in Z-axis)

Post Processing

None


Image: K. Small

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000053
Stage	Western blot
Instrument	Azure Sapphire FL
Microscope configuration	Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)
Structured source metadata values	11
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	B4DF57E56A308FEFC293951F62082C04171A7FEFC54F21EBE0AEE02590300896
Method PDF SHA-256	BAFCC67EB85861EF7000DBCC56BB3E980A488F194C39845DDC8F53E8D2113F4D

ORIGINAL IMAGE EVIDENCE / WIDEFIELD

A

B

	Channel 1	Channel 2	Channel 3
Reflector	AF532-49014	50 Cy 5	49 DAPI
Beam Splitter	550	660	395
Filter Ex. Wavelength	515 - 545	625 - 655	335 - 383
Filter Em. Wavelength	555 - 595	665 - 715	420 - 470
Contrast Method	Fluorescence	Fluorescence	Fluorescence
Light Source	X-Cite Xylis Lamp	X-Cite Xylis Lamp	X-Cite Xylis Lamp
Light Source Intensity	25.00 %	25.00 %	10.00 %
Channel Name	AF532	AF647	DAPI
Channel Description			
Dye Name	AF532	AF647	DAPI
Channel Color			
Excitation Wavelength	528	653	353
Emission Wavelength	553	668	465
Effective NA	1.4	1.4	1.4
Imaging Device	Axiocam 705	Axiocam 705	Axiocam 705
Camera Adapter	1x Camera Adapter	1x Camera Adapter	1x Camera Adapter
Exposure Time	250 ms	250 ms	25 ms
Depth of Focus	0.86 μm	1.03 μm	0.72 μm
Binning Mode	2,2	2,2	2,2

Widefield

Baseline widefield imaging of the PARP1 Alexa Fluor™ 532 Direct Conjugate, generated using the Rabbit anti-Human PARP1 Recombinant Monoclonal antibody, acquired following 2 μM etoposide treatment, demonstrates robust nuclear localization and DNA damage-associated signal. These features are equivalently observed in the unconjugated PARP1 reference condition (Cat. ID# PA-RR-000003), with Zeiss ZEN image acquisition settings displayed for transparency and held constant across conjugation states. Collectively, these data meet CDx go/no-go criteria for baseline equivalency and support downstream analytical and enrichment applications within Identifyn®'s step-down imaging methodology.

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn™ SRM Alexa Fluor™ 532, IgG - BLR018M
Cat ID# - DC 000053

EVIDENCE TIER	Widefield
INSTRUMENT	ZEISS Axio Imager.M1
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 532; 555; 647
METADATA VALUES	37

STRUCTURED ACQUISITION METADATA / WIDEFIELD

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705 mono, and X-Cite Xylis Lamp, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63x/1.40 Oil M27
Channels	3
Scaling (Per Pixel)	0.110µm X 0.110µm
Image size (Pixels)	725 X 644
Image Size (Scaled)	79.40 µm X 644 µm
Bit Depth	14 Bit
Image Center Position	X: 0.00µm, Y: 0.00µm
ROI Center Offset	X: 0.00µm, Y: 0.00µm
Reflector	AF532-49014 50 Cy 5 49 DAPI
Beam Splitter	550 660 395
Filter Excitation Wavelength	515 - 545 625 - 655 335 - 383
Filter Emission Wavelength	555 - 595 665 - 715 420 - 470
Contrast Method	Fluorescence
Light source	X-Cite Xylis Lamp Intensity @25.00% X-Cite Xylis Lamp Intensity @10.00%
Exposure Time	250 ms 25 ms
Excitation Wavelength	528 653 353
Emission Wavelength	553 668 465
Effective NA	1.4
Imaging Device	Axiocam 705
Camera Adapter	1X
Depth of Focus	0.86 µm 1.03 µm 0.72 µm
Binning Mode	2,2

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L)
Cat ID# - SA 000059

PARP1 (Poly(ADP-ribose) polymerase 1) has a multifaceted function in the detection of DNA damage, including single-strand breaks (SSBs) and double-strand breaks (DSBs) via binding at the site of damage. PARP1 catalyzes the synthesis of poly(ADP-ribose) (PAR) chains using NAD⁺; the resulting PAR chains signal the recruitment of DNA repair proteins to the break, initiating the repair process. PARP1, as well as PARP2, are retained in nucleoli via interaction with the multifunctional nucleolar hub protein nucleophosmin, which is implicated in multiple steps of ribosome biogenesis, including rDNA transcription and elongation, as well as rRNA processing.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS.

The second primary antibody, Nuclear Pore Complex Monoclonal Antibody (39C7), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705 mono, and X-Cite Xylis Lamp, was used with the following parameters for each dye:

Single Plane (2D) - Split View Panel A

Channels = 3

Scaling (Per Pixel) = 0.110µm X 0.110µm

Image size (Pixels) = 725 X 644

Image Size (Scaled) = 79.40 µm X 644 µm

Bit Depth = 14 Bit

Image Center Position = X: 0.00µm, Y: 0.00µm

ROI Center Offset = X: 0.00µm, Y: 0.00µm

532-

Reflector = AF532-49014
 Beam Splitter = 550
 Filter Excitation Wavelength = 515 - 545
 Filter Emission Wavelength = 555 - 595
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @25.00%
 Exposure Time = 250 ms
 Excitation Wavelength = 528
 Emission Wavelength = 553
 Effective NA = 1.4
 Imaging Device = Axiocam 705
 Camera Adapter = 1X
 Depth of Focus = 0.86 μm
 Binning Mode = 2,2

647 -

Reflector = 50 Cy 5
 Beam Splitter = 660
 Filter Excitation Wavelength = 625 - 655
 Filter Emission Wavelength = 665 - 715
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @25.00%
 Exposure Time = 250 ms
 Excitation Wavelength = 653
 Emission Wavelength = 668
 Effective NA = 1.4
 Imaging Device = Axiocam 705
 Camera Adapter = 1X
 Depth of Focus = 1.03 μm
 Binning Mode = 2,2

DAPI -

Reflector = 49 DAPI
 Beam Splitter = 395
 Filter Excitation Wavelength = 335 - 383
 Filter Emission Wavelength = 420 - 470
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @10.00%
 Exposure Time = 25 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Imaging Device = Axiocam 705
 Camera Adapter = 1X
 Depth of Focus = 0.72 μm
 Binning Mode = 2,2

Split View - Panel A

The top-left panel shows Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate. The top-right panel shows monoclonal nuclear pore complex (NPC), visualized using Identifyn® SRM Alexa Fluor™ 647 affinity-purified donkey anti-mouse IgG (H+L). The lower-left panel shows DAPI staining, and the bottom-right panel shows the three-channel merged image.

Zen Acquisition Info Tab - Panel B

Shown is the acquisition of the 532 nm channel (H2AX Alexa Fluor™ 532 Direct Conjugate) using an X-Cite Xylis illumination source operated at 25% output and detected with an Axiocam 705 camera using a 200 ms exposure time. Additionally shown is the 647 channel (NPC/647) using 25% input with a 200ms exposure time.

Summary

Baseline widefield imaging of the PARP1 Alexa Fluor™ 532 Direct Conjugate, generated using the Rabbit anti-Human PARP1 Recombinant Monoclonal antibody, acquired following 2 µM etoposide treatment, demonstrates robust nuclear localization and DNA damage-associated signal. These features are equivalently observed in the unconjugated PARP1 reference condition (Cat. ID# PA-RR-000003), with Zeiss ZEN image acquisition settings displayed for transparency and held constant across conjugation states. Collectively, these data meet CDx go/no-go criteria for baseline equivalency and support downstream analytical and enrichment applications within Identifyn®'s step-down imaging methodology.

Image Corrections

NONE

Image by E.F.Simon


Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000053
Stage	Widefield
Instrument	ZEISS Axio Imager.M1
Microscope configuration	Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705 mono, and X-Cite Xylis Lamp, was used with the following parameters for each dye:
Structured source metadata values	37
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	D347ABC0E519C78DC6E60DC3B66F46CFF2FF0E609D35D949984F2C2C1ABEC56D
Method PDF SHA-256	7DCC9719F0B2F461A60FCBC68A523FBA8425FA8B9CD387B8452E46197EFC05F4

ORIGINAL IMAGE EVIDENCE / WIDEFIELD SIM - APOTOME



EVIDENCE TIER	Widefield SIM - Apotome
INSTRUMENT	ZEISS Axio Observer.Z1/7 Apotome
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 532; 555; 647
METADATA VALUES	56

STRUCTURED ACQUISITION METADATA / WIDEFIELD SIM - APOTOME

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:
Objective	EC Plan-Neofluar 63X/1.25 oil M27
Channels	3
Scaling (Per Pixel)	0.143 μ m X 0.143 μ m X 0.100 μ m
Image size (Pixels)	535 X 417
Image Size (Scaled)	76.43 μ m X 59.57 μ m
Bit Depth	14 Bit
Image Center Position	X: 63.21 μ m, Y: 46.97 μ m
Stage Position	X: 63.21 μ m, Y: 46.97 μ m
ROI Center Offset	X: 1.14 μ m, Y: 0.00 μ m
Z-Stack	98 Slices (9.7 μ m)
Reflector	AF532-49014 50 Cy 5 49 DAPI
Beam Splitter	550 660 395
Filter Excitation Wavelength	515 - 545 625 - 655 335 - 383
Filter Emission Wavelength	555 - 595 665 - 715 420 - 470
Contrast Method	Fluorescence
Light source	X-Cite Xylis Lamp Intensity @20.00% X-Cite Xylis Lamp Intensity @10.00%
Exposure Time	200 ms 25 ms
Excitation Wavelength	528 653 353
Emission Wavelength	553 668 465
Effective NA	1.4 1.25
Imaging Device	Axiocam 807 Mono
Camera Adapter	1X
Depth of Focus	1.07 μ m 1.30 μ m 0.90 μ m
Binning Mode	2,2
Apotome SIM Image Processing	Display Mod" "Optical Sectioning" - enabled correction selected.
3D Maximum Projection	Precise
Appearance	Threshold 13% all channels, Ramp 0% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 100%, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
X-Y	Activate was selected, textured opaque was chosen, and a red outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

FIELD	SOURCE-DOCUMENTED VALUE
Position of Clip	11% 13% -17%
Clip X	0°
Clip Z	0°
X-Z	Activate was selected, textured opaque was chosen, and a yellow outline of the clipped perimeter was shown; Clip Front was selected, and Clip Back was NOT selected.
Clip Y	-45° 45°
Y-Z	Activate was selected, textured opaque was chosen, and a blue outline of the clipped perimeter was shown; Clip Front was selected, and Clip Back was NOT selected.

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L)
Cat ID# - SA 000059

PARP1 (Poly(ADP-ribose) polymerase 1) has a multifaceted function in the detection of DNA damage, including single-strand breaks (SSBs) and double-strand breaks (DSBs) via binding at the site of damage. PARP1 catalyzes the synthesis of poly(ADP-ribose) (PAR) chains using NAD⁺; the resulting PAR chains signal the recruitment of DNA repair proteins to the break, initiating the repair process. PARP1, as well as PARP2, are retained in nucleoli via interaction with the multifunctional nucleolar hub protein nucleophosmin, which is implicated in multiple steps of ribosome biogenesis, including rDNA transcription and elongation, as well as rRNA processing.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS.

The second primary antibody, Nuclear Pore Complex Monoclonal Antibody (39C7), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:

Single Plane (2D) - Panel A, Z-plane 47 of 98

Channels = 3

Scaling (Per Pixel) = 0.143 µm X 0.143 µm X 0.100 µm

Image size (Pixels) = 535 X 417

Image Size (Scaled) = 76.43 µm X 59.57 µm

Bit Depth = 14 Bit

Image Center Position = X: 63.21 µm, Y: 46.97 µm

Stage Position = X: 63.21 µm, Y: 46.97 µm

ROI Center Offset = X: 1.14 µm, Y: 0.00 µm

Z Stack - (3D) - Panel B

Z-Stack = 98 Slices (9.7 µm)

Channels = 3

Scaling (Per Pixel) = 0.143 µm X 0.143 µm X 0.100 µm

Image size (Pixels) = 535 X 417

Image Size (Scaled) = 76.43 µm X 59.57 µm

Bit Depth = 14 Bit

Image Center Position = X: 63.21 µm, Y: 46.97 µm

Stage Position = X: 63.21 µm, Y: 46.97 µm

ROI Center Offset = X: 1.14 µm, Y: 0.00 µm

532 -

Reflector = AF532-49014

Beam Splitter = 550

Filter Excitation Wavelength = 515 - 545

Filter Emission Wavelength = 555 - 595

Contrast Method = Fluorescence

Light source = X-Cite Xylis Lamp Intensity @20.00%

Exposure Time = 200 ms

Excitation Wavelength = 528

Emission Wavelength = 553

Effective NA = 1.4

Imaging Device = AxioCam 807 Mono

Camera Adapter = 1X

Depth of Focus = 1.07 µm

Binning Mode = 2,2

647 -

Reflector = 50 Cy 5

Beam Splitter = 660

Filter Excitation Wavelength = 625 - 655

Filter Emission Wavelength = 665 - 715

Contrast Method = Fluorescence

Light source = X-Cite Xylis Lamp Intensity @20.00%

Exposure Time = 200 ms

Excitation Wavelength = 653

Emission Wavelength = 668

Effective NA = 1.25

Imaging Device = AxioCam 807 Mono

Camera Adapter = 1X

Depth of Focus = 1.30 µm

Binning Mode = 2,2

DAPI - (not shown)

Reflector = 49 DAPI
 Beam Splitter = 395
 Filter Excitation Wavelength = 335 - 383
 Filter Emission Wavelength = 420 - 470
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @10.00%
 Exposure Time = 25 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Imaging Device = AxioCam 807 Mono
 Camera Adapter = 1X
 Depth of Focus = 0.90 μ m
 Binning Mode = 2,2

Post Processing - SIM Processing

Apotome SIM Image Processing = Display Mod" "Optical Sectioning" - enabled correction selected.
 Normalization, Phase Corrector, Fourier Filter, and Deconvolution were NOT used.

2-Dimensional View - Panel A

Shows Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M) and Mouse monoclonal nuclear pore complex primary with Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L). Taken from Z-plane 47 of 98.

3D Image Processing - Panel B

3D Maximum Projection = Precise
 Appearance = Threshold 13% all channels, Ramp 0% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 100%, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Clipping Image Process - Panel B

Clipping planes are introduced in Panel B to highlight the PARP1/532 direct conjugate within the nuclear compartment, using the NPC counterstained with 647 as an exclusionary marker. This is achieved by adjusting or cutting away portions of the image using a triple-wedge cut, removing parts of the image in the X-Y, X-Z, and Y-Z planes.

X-Y = Activate was selected, textured opaque was chosen, and a red outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = 11%

Clip X = 0°

Clip Z = 0°

X-Z = Activate was selected, textured opaque was chosen, and a yellow outline of the clipped perimeter was shown; Clip Front was selected, and Clip Back was NOT selected.

Position of Clip = 13%

Clip Y = -45°

Clip Z = 0°

Y-Z = Activate was selected, textured opaque was chosen, and a blue outline of the clipped perimeter was shown; Clip Front was selected, and Clip Back was NOT selected.

Position of Clip = -17%

Clip Y = 45°

Clip Z = 0°

Panel B shows a wedge-cut view displayed from the top with a slight backward rotation, highlighting the region of wedge removal and revealing the full depth of signal across the z-stack. This orientation provides an alternative view of wedge geometry and confirms the intranuclear localization of PARP1/532, as indicated by its colocalization relative to the NPC counterstained with the Identifyn® SRM Alexa Fluor™ 647.

Summary

Widefield SIM (Apotome 3D) imaging shows a 2D optical section (Panel A) of PARP1 Alexa Fluor™ 532 (yellow), detected using the Rabbit anti-Human PARP1 Recombinant Monoclonal antibody directly conjugated with Identifyn® SRM Alexa Fluor™ 532, together with the Nuclear Pore Complex (NPC) labeled using a Mouse monoclonal NPC antibody detected with Identifyn® SRM Alexa Fluor™ donkey anti-mouse 647 (red). A corresponding three-dimensional maximum-intensity projection (Panel B) demonstrates compartment-restricted nuclear localization of PARP1 following etoposide-induced DNA damage, with the NPC signal serving as an exclusionary boundary marker to confirm nuclear inclusion. Improved optical sectioning relative to standard widefield imaging is observed without asserting sub-nuclear structural resolution. PARP1 signal distribution and intensity are consistent with the unconjugated PARP1 reference condition (Cat. ID# PA-RR-000003), satisfying CDx go/no-go criteria for conjugation equivalency, localization specificity, and analytical suitability within Identifyn®'s step-down imaging methodology.

Image Corrections

NONE

Image by E.F.Simon

Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.


SOURCE PROVENANCE

Catalog	DC-000053
Stage	Widefield SIM - Apotome
Instrument	ZEISS Axio Observer.Z1/7 Apotome
Microscope configuration	Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:
Structured source metadata values	56
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT

SOURCE PROVENANCE

Cell model	HeLa
Image SHA-256	CE77DBCE5EE34531408770D51C31E21B4AC7C595596AD9701A35DBC3E7CE958B
Method PDF SHA-256	0B734065CAD892077DDD3BF86DBC65C5F29AC82739DBA9391E799701442C79B8

ORIGINAL IMAGE EVIDENCE / CONFOCAL



EVIDENCE TIER	Confocal
INSTRUMENT	ZEISS LSM 980
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 532; 647
METADATA VALUES	59

STRUCTURED ACQUISITION METADATA / CONFOCAL

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	121 Slices (3.6 μm)
Channels	3
Scaling (Per Pixel)	0.059 μm X 0.059 μm X 0.030 μm
Image size (Pixels)	1121 X 1121
Image Size (Scaled)	65.93 μm X 65.93 μm
Bit Depth	16 Bit
Image Center Position	X: 92.14 μm , Y: 54.19 μm
Stage Position	X: 92.14 μm , Y: 54.19 μm
ROI Center Offset	X: 0.00 μm , Y: 0.00 μm
Reflector	None
Contrast Method	Fluorescence
Pinhole	1.00AU / 53 μm 1.00AU / 66 μm 1.00AU / 43 μm
LMS MBS	MBS 488/639 & 405/730
LMS SBS	Mirror SBS SP 615
Laser Wavelength	488nm @0.8% 639nm @0.10% 405nm @0.2%
Laser Blanking	Enabled
Scan Mode	Frame
Scan Zoom	2.0
Rotation	0°
Sampling	1.20
Pixel Time	0.30 μs
Frame Time	3.69s
LSM Scan Speed	12
Scan Direction	Bidirectional
Line Step	1
Averaging	8
Excitation Wavelength	528 653 353
Emission Wavelength	553 668 465
Effective NA	1.4
Detection Wavelength	530 - 580 643 - 740 440 - 470
Imaging Device	GaAsP-Pmt3 GaAsP-Pmt3

FIELD	SOURCE-DOCUMENTED VALUE
Detector Type	GaAsP-PMT GaAsP-PMT
Detector Gain	550 V 500 V
Detector Offset	0
Detector Digital Gain	1.0
LSM Plus Mode	Filter: 6.8 (3D, Auto) Filter: 7.5 (3D, Auto) Filter: 6.3 (3D, Auto)
Profile Definition	Stroke Thickness 1, Normal View, Grid Distance 1
Profile View	X and Y axes were set to Auto - X (Min 0.0 and Max 67.4), Y (Min 0 and Max 5600)
3D Maximum Projection	Precise
Appearance	Threshold 1% all channels, Ramp 0% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 100%, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L)
Cat ID# - SA 000059

PARP1 (Poly(ADP-ribose) polymerase 1) has a multifaceted function in the detection of DNA damage, including single-strand breaks (SSBs) and double-strand breaks (DSBs) via binding at the site of damage. PARP1 catalyzes the synthesis of poly(ADP-ribose) (PAR) chains using NAD⁺; the resulting PAR chains signal the recruitment of DNA repair proteins to the break, initiating the repair process. PARP1, as well as PARP2, are retained in nucleoli via interaction with the multifunctional nucleolar hub protein nucleophosmin, which is implicated in multiple steps of ribosome biogenesis, including rDNA transcription and elongation, as well as rRNA processing.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS.

The second primary antibody, Nuclear Pore Complex Monoclonal Antibody (39C7), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:

Single Plane (2D) - Panel A, Z-plane 32 of 121

Z Stack - (3D) - Panel D

Z-Stack = 121 Slices (3.6 µm)

Channels = 3

Scaling (Per Pixel) = 0.059 µm X 0.059 µm X 0.030 µm

Image size (Pixels) = 1121 X 1121

Image Size (Scaled) = 65.93 µm X 65.93 µm

Bit Depth = 16 Bit

Image Center Position = X: 92.14 µm, Y: 54.19 µm

Stage Position = X: 92.14 µm, Y: 54.19 µm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

532 -

Reflector = None

Contrast Method = Fluorescence

Pinhole = 1.00AU / 53µm

LMS MBS = MBS 488/639 & 405/730

LMS SBS = Mirror

Laser Wavelength = 488nm @0.8%

Laser Blanking = Enabled

Scan Mode = Frame

Scan Zoom = 2.0

Rotation = 0°

Sampling = 1.20

Pixel Time = 0.30µs

Frame Time = 3.69s

LSM Scan Speed = 12

Scan Direction = Bidirectional

Line Step = 1

Averaging = 8

Excitation Wavelength = 528

Emission Wavelength = 553

Effective NA = 1.4

Detection Wavelength = 530 - 580

Imaging Device = GaAsP-Pmt3

Detector Type = GaAsP-PMT

Detector Gain = 550 V

Detector Offset = 0

Detector Digital Gain = 1.0

LSM Plus Mode = Filter: 6.8 (3D, Auto)

647 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 1.00AU / 66µm
 LMS MBS = MBS 488/639 & 405/730
 LMS SBS = SBS SP 615
 Laser Wavelength = 639nm @0.10%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.0
 Rotation = 0°
 Sampling = 1.20
 Pixel Time = 0.30µs
 Frame Time = 3.69s
 LSM Scan Speed = 12
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 8
 Excitation Wavelength = 653
 Emission Wavelength = 668
 Effective NA = 1.4
 Detection Wavelength = 643 - 740
 Imaging Device = GaAsP-Pmt3
 Detector Type = GaAsP-PMT
 Detector Gain = 550 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 LSM Plus Mode = Filter: 7.5 (3D, Auto)

DAPI - (Not Shown in Data Analysis)

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 1.00AU / 43µm
 LMS MBS = MBS 488/639 & 405/730
 LMS SBS = Mirror
 Laser Wavelength = 405nm @0.2%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.0
 Rotation = 0°
 Sampling = 1.20
 Pixel Time = 0.30µs
 Frame Time = 3.69s
 LSM Scan Speed = 12
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 8
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Detection Wavelength = 440 - 470
 Imaging Device = GaAsp-Pmt3
 Detector Type = GaAsp-PMT
 Detector Gain = 500 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 LSM Plus Mode = Filter: 6.3 (3D, Auto)

Profile View Display - Panels A and B

Profile Definition = Stroke Thickness 1, Normal View, Grid Distance 1
 Profile View = X and Y axes were set to Auto - X (Min 0.0 and Max 67.4), Y (Min 0 and Max 5600)

This panel presents a line profile measured along an arbitrary axis at z-plane 32, quantifying the fluorescence intensity of the PARP1 Alexa Fluor™ 532 direct conjugate, which reaches a maximum of 5,333 A.U. A corresponding line profile of the NPC/647 signal delineates the nuclear envelope, confirming that the PARP1 532 signal is confined within the nuclear compartment. In comparison, the PARP1 532 direct conjugate shows higher brightness than NPC staining.

Zen Acquisition Info Tab - Panel C

Shown is the acquisition of the 532 channel (PARP1 532 direct conjugate) using the 488nm laser at 0.8%, demonstrating that the signal is not artificially produced by image acquisition settings.

3D Image Processing

3D Maximum Projection = Precise

Appearance = Threshold 1% all channels, Ramp 0% all channels, Maximum 100% all channels

Background = Black

Light = Brightness 100%, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Confocal imaging of PARP1 using the Identifyn™ SRM Alexa Fluor™ 532 Direct Conjugate, counterstained with a Nuclear Pore Complex (NPC) marker detected in the 647 nm channel as an exclusionary boundary control, demonstrates damage-dependent nuclear recruitment following etoposide treatment. A representative two-dimensional optical section acquired at z-plane 32 reveals a robust intranuclear PARP1 signal with a maximum fluorescence intensity of 5,333.00 A.U., consistent with pharmacologically induced DNA damage response activation.

Acquisition transparency is provided via Zeiss ZEN metadata, confirming low-power imaging conditions (488 nm excitation at 0.8%), ensuring that observed signal intensity reflects a genuine biological response rather than acquisition-driven amplification. A corresponding three-dimensional maximum-intensity projection with orthogonal views (X-Y, X-Z, Y-Z) further confirms confinement of PARP1 signal within the nuclear compartment relative to the NPC boundary, satisfying CCR enrichment criteria for functional recruitment.

At confocal resolution, PARP1 is observed as punctate and aggregated intranuclear signal consistent with PARylation-associated accumulation; however, these features remain structurally unresolved and likely represent resolution-limited composites. Definitive characterization of PARP1 sub-assemblies and chromatin-associated organization requires progression beyond the diffraction limit into super-resolution modalities (Airyscan, SIM, SIM², and STORM), where structural organization can be adjudicated within the CCR step-down framework.

Collectively, these confocal data satisfy CDx enrichment and go/no-go criteria by confirming damage-dependent PARP1 recruitment under low-power conditions, while appropriately deferring mechanistic structural interpretation to higher-resolution gates.

Image Corrections

NONE

Image by J.K. Bennett

Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.


SOURCE PROVENANCE

Catalog	DC-000053
Stage	Confocal
Instrument	ZEISS LSM 980
Microscope configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:

SOURCE PROVENANCE

Structured source metadata values	59
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	02C659102B981DA3B4CF815C26EE2C90EDEC20E7D3D008B9FEE2D2EE6EB920AB
Method PDF SHA-256	B8F71611DCE3AA1A688391230F5E50EF283FF6BFEB1DEFCE6DA2B0D14EBE14B8

ORIGINAL IMAGE EVIDENCE / AIRYSCAN



EVIDENCE TIER	Airyscan
INSTRUMENT	ZEISS LSM 980
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 532; 647; 750
METADATA VALUES	59

STRUCTURED ACQUISITION METADATA / AIRYSCAN

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Channels	3
Scaling (Per Pixel)	0.03530 μm X 0.03530 μm X 0.03000 μm
Image size (Pixels)	607 X 798 241 X 178
Image Size (Scaled)	21.42 μm X 28.14 μm 8.51 μm X 6.28 μm
Bit Depth	16 Bit
Image Center Position	X: 92.50 μm , Y: 53.72 μm
Stage Position	X: 92.50 μm , Y: 53.72 μm
ROI Center Offset	X: 0.00 μm , Y: 0.00 μm
Z-Stack	141 Slices (4.2 μm)
Reflector	None
Contrast Method	Fluorescence
Pinhole	5.00AU / 265 μm 5.00AU / 328 μm 5.00AU / 214 μm
LMS MBS	MBS 488/639 & 405/730
LMS SBS	SBS LP 525 SBS LP 640 SBS SP 505
Laser Wavelength	488nm @1.20% 639nm @0.10% 405nm @0.2%
Laser Blanking	Enabled
Scan Mode	Frame
Scan Zoom	2.4
Rotation	0°
Sampling	2.00
Pixel Time	0.66 μs
Frame Time	7.82s
LSM Scan Speed	7
Scan Direction	Bidirectional
Line Step	1
Averaging	4
Excitation Wavelength	528 653 353
Emission Wavelength	553 668 465
Effective NA	1.4
Detection Wavelength	530 - 582 643 - 735 422 - 477
Imaging Device	Airyscan

FIELD	SOURCE-DOCUMENTED VALUE
Detector Type	GaAsP-PMT GaAsP-PMT
Detector Gain	800 V 750 V
Detector Offset	0
Detector Digital Gain	1.0
Super Resolution	9.2 (3D, Auto) 9.8 (3D, Auto) 8.1 (3D, Auto)
3D Volume Projection	Precise
Appearance	Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L)
Cat ID# - SA 000059

PARP1 (Poly(ADP-ribose) polymerase 1) has a multifaceted function in the detection of DNA damage, including single-strand breaks (SSBs) and double-strand breaks (DSBs) via binding at the site of damage. PARP1 catalyzes the synthesis of poly(ADP-ribose) (PAR) chains using NAD⁺; the resulting PAR chains signal the recruitment of DNA repair proteins to the break, initiating the repair process. PARP1, as well as PARP2, are retained in nucleoli via interaction with the multifunctional nucleolar hub protein nucleophosmin, which is implicated in multiple steps of ribosome biogenesis, including rDNA transcription and elongation, as well as rRNA processing.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS.

The second primary antibody, Nuclear Pore Complex Monoclonal Antibody (39C7), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Nuclear Pore Complex(NPC), Identifyn® SRM Alexa Fluor™ 647 Staining not represented in data analysis

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:

Single Plane (2D) - Panel A, Z-plane 71 of 141

Channels = 3

Scaling (Per Pixel) = 35.30 nm X 35.30 nm X 30.00 nm

Image size (Pixels) = 607 X 798

Image Size (Scaled) = 21.42 µm X 28.14 µm

Bit Depth = 16 Bit

Image Center Position = X: 92.50 µm, Y: 53.72 µm

Stage Position = X: 92.50 µm, Y: 53.72 µm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D) - Panels B, C, and D

Z-Stack = 141 Slices (4.2 µm)

Channels = 3

Scaling (Per Pixel) = 35.30 nm X 35.30 nm X 30.00 nm

Image size (Pixels) = 241 X 178

Image Size (Scaled) = 8.51 µm X 6.28 µm

Bit Depth = 16 Bit

Image Center Position = X: 92.50 µm, Y: 53.72 µm

Stage Position = X: 92.50 µm, Y: 53.72 µm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

532 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 5.00AU / 265µm
 LMS MBS = MBS 488/639 & 405/730
 LMS SBS = SBS LP 525
 Laser Wavelength = 488nm @1.20%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.4
 Rotation = 0°
 Sampling = 2.00
 Pixel Time = 0.66µs
 Frame Time = 7.82s
 LSM Scan Speed = 7
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 4
 Excitation Wavelength = 528
 Emission Wavelength = 553
 Effective NA = 1.4
 Detection Wavelength = 530 - 582
 Imaging Device = Airyscan
 Detector Type = GaAsP-PMT
 Detector Gain = 800 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 Super Resolution = 9.2 (3D, Auto)

647 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 5.00AU / 328µm
 LMS MBS = MBS 488/639 & 405/730
 LMS SBS = SBS LP 640
 Laser Wavelength = 639nm @0.10%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.4
 Rotation = 0°
 Sampling = 2.00
 Pixel Time = 0.66µs
 Frame Time = 7.82s
 LSM Scan Speed = 7
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 4
 Excitation Wavelength = 653
 Emission Wavelength = 668
 Effective NA = 1.4
 Detection Wavelength = 643 - 735
 Imaging Device = Airyscan
 Detector Type = GaAsP-PMT
 Detector Gain = 800 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 Super Resolution = 9.8 (3D, Auto)

DAPI - (Not Shown in Data Analysis)

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 5.00AU / 214µm
 LMS MBS = MBS 488/639 & 405/730
 LMS SBS = SBS SP 505
 Laser Wavelength = 405nm @0.2%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.4
 Rotation = 0°
 Sampling = 2.00
 Pixel Time = 0.66µs
 Frame Time = 7.82s
 LSM Scan Speed = 7
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 4
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Detection Wavelength = 422 - 477
 Imaging Device = Airyscan
 Detector Type = GaAsP-PMT
 Detector Gain = 750 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 Super Resolution = 8.1 (3D, Auto)

2-Dimensional View - Panel A

Panel A - Shows Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M) at Z-plane 71 of 141 in default yellow with a white box signifying the region of interest to be further examined in panels B, C, and D.

3D Image Processing - Panels B, C, and D

3D Volume Projection = Precise
 Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Airyscan imaging of the PARP1 Alexa Fluor™ 532 Direct Conjugate following etoposide-induced DNA damage is shown as a 2D optical section at Z-plane 71 of 141 (Panel A), with a region of interest (ROI, white box) selected for downstream structural evaluation. Panel B presents a 3D volume projection of DAPI alone, delineating chromatin architecture and nuclear exclusion zones, including the nucleolus. Panel C shows the corresponding PARP1 volume projection, revealing damage-associated PARP1 aggregates resolved at Airyscan resolution. Panel D displays the merged PARP1/DAPI volume, demonstrating consistent localization of PARP1 aggregates within and surrounding nucleolar chromatin boundaries. At this resolution (~120 nm lateral), PARP1 aggregates likely represent unresolved composites associated with active PARYlation states; definitive substructure requires progression to higher-resolution modalities (SIM/SIM², STORM). Collectively, these data satisfy CCR enrichment and transitional-resolution CDx criteria, confirming damage-dependent nuclear recruitment while appropriately constraining mechanistic interpretation pending higher-resolution validation.

Image Corrections

NONE

Image by J.K. Bennett


Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000053
Stage	Airyscan
Instrument	ZEISS LSM 980
Microscope configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Structured source metadata values	59
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	99059AC731C4365B76A99326096F4A9A2E6328B6BBBA546D6A68F238F785339D
Method PDF SHA-256	02D7DE4383B706F787B7C1BFDF1795A0D40B4E0F1B9E1A926ACCCD2D93337222

ORIGINAL IMAGE EVIDENCE / SIM



EVIDENCE TIER	SIM
INSTRUMENT	ZEISS Elyra
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 532; 647
METADATA VALUES	65

STRUCTURED ACQUISITION METADATA / SIM

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	142 Slices (4.23 μm)
Channels	3
Scaling (Per Pixel)	0.02111 μm X 0.02111 μm X 0.03000 μm
Image size (Pixels)	416 X 462 184 X 164
Image Size (Scaled)	8.78 μm X 9.75 μm 3.88 μm X 3.46 μm
Bit Depth	16 Bit
Image Center Position	X: 58.05 μm , Y: 35.69 μm
Stage Position	X: 58.05 μm , Y: 35.69 μm
ROI Center Offset	X: 0.00 μm , Y: 0.00 μm
Reflectors	SR HE LBF 405/488/561/642
Beam Splitter	405, 488, 561, 642
Filter Excitation Wavelength	200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method	Fluorescence
Illumination Wavelength	488 640 405
Channel Name	T2-Axiocam 820 #2-SR T1-Axiocam 820 -SR T3-Axiocam 820 #2-SR
Excitation Wavelength	488 640 405
Emission Wavelength	525 729 445
Effective NA	1.4
Detection Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800 503 - 546, 657 - 800
Imaging Device	Axiocam 820 #2 Axiocam 820
Camera Adapter	1X
Exposure Time	150 ms 100 ms 50 ms
Depth of Focus	0.81 μm 1.13 μm 0.69 μm
Binning Mode	2,2
Laser Wavelength	488 nm @4.0% 640 nm @2.0% 405 nm @2.0%
Frame Time	1.59 s 1.33 s
Scan Direction	Unidirectional
Grating	36.5 μm grating 27.5 μm grating
Processing Dimension	3D
Result Sampling	4

FIELD	SOURCE-DOCUMENTED VALUE
Process Sampling	N/A
Algorithm	SIM
SIM Settings	Standard
Auto Sharpness	On
Sharpness	9.8786 (647), 11.1960 (532), 9.9820 (DAPI)
Fitting	Fast Fit
Normalization	Auto
3D Volume Projection	Precise
Appearance	Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 125%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected) View Angle 35°, Scale Z 35%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L)
Cat ID# - SA 000059

PARP1 (Poly(ADP-ribose) polymerase 1) has a multifaceted function in the detection of DNA damage, including single-strand breaks (SSBs) and double-strand breaks (DSBs) via binding at the site of damage. PARP1 catalyzes the synthesis of poly(ADP-ribose) (PAR) chains using NAD⁺; the resulting PAR chains signal the recruitment of DNA repair proteins to the break, initiating the repair process. PARP1, as well as PARP2, are retained in nucleoli via interaction with the multifunctional nucleolar hub protein nucleophosmin, which is implicated in multiple steps of ribosome biogenesis, including rDNA transcription and elongation, as well as rRNA processing.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS.

The second primary antibody, Nuclear Pore Complex Monoclonal Antibody (39C7), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Nuclear Pore Complex(NPC), Identifyn® SRM Alexa Fluor™ 647 Staining not represented in data analysis

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:

Z Stack - (3D) - Panel A, and Panel B, Z-plane 56 of 142

Z-Stack = 142 Slices (4.23 µm)

Channels = 3

Scaling (Per Pixel) = 21.11 nm X 21.11 nm X 30.00 nm

Image size (Pixels) = 416 X 462

Image Size (Scaled) = 8.78 µm X 9.75 µm

Bit Depth = 16 Bit

Image Center Position = X: 58.05 µm, Y: 35.69 µm

Stage Position = X: 58.05 µm, Y: 35.69 µm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D) - Panel C

Z-Stack = 142 Slices (4.23 µm)

Channels = 3

Scaling (Per Pixel) = 21.11 nm X 21.11 nm X 30.00 nm

Image size (Pixels) = 184 X 164

Image Size (Scaled) = 3.88 µm X 3.46 µm

Bit Depth = 16 Bit

Image Center Position = X: 58.05 µm, Y: 35.69 µm

Stage Position = X: 58.05 µm, Y: 35.69 µm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

532 -

Reflectors = SR HE LBF 405/488/561/642

Beam Splitter = 405, 488, 561, 642

Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642

Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800

Contrast Method = Fluorescence

Illumination Wavelength = 488

Channel Name = T2-Axiocam 820 #2-SR

Excitation Wavelength = 488

Emission Wavelength = 525

Effective NA = 1.4

Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800

Imaging Device = Axiocam 820 #2

Camera Adapter = 1X

Exposure Time = 150 ms

Depth of Focus = 0.81 µm

Binning Mode = 2,2

Laser Wavelength = 488 nm @4.0%

Frame Time = 1.59 s

Scan Direction = Unidirectional

Grating = 36.5 µm grating

647 - (Not Shown in Data Analysis)

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 640
 Channel Name = T1-Axiocam 820 -SR
 Excitation Wavelength = 640
 Emission Wavelength = 729
 Effective NA = 1.4
 Detection Wavelength = 503 - 546, 657 - 800
 Imaging Device = Axiocam 820
 Camera Adapter = 1X
 Exposure Time = 100 ms
 Depth of Focus = 1.13 μm
 Binning Mode = 2,2
 Laser Wavelength = 640 nm @2.0%
 Frame Time = 1.33 s
 Scan Direction = Unidirectional
 Grating = 27.5 μm grating

DAPI -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 405
 Channel Name = T3-Axiocam 820 #2-SR
 Excitation Wavelength = 405
 Emission Wavelength = 445
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820 #2
 Camera Adapter = 1X
 Exposure Time = 50 ms
 Depth of Focus = 0.69 μm
 Binning Mode = 2,2
 Laser Wavelength = 405 nm @2.0%
 Frame Time = 1.33 s
 Scan Direction = Unidirectional
 Grating = 36.5 μm grating

Post Processing - SIM Processing

Processing Dimension = 3D

Result Sampling = 4

Process Sampling = N/A

Channels = 3

Algorithm = SIM

SIM Settings = Standard

Auto Sharpness = On

Sharpness = 9.8786 (647), 11.1960 (532), 9.9820 (DAPI)

Fitting = Fast Fit

Normalization = Auto

3D Image Processing - Panels A

3D Volume Projection = Precise

Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels

Background = Black

Light = Brightness 125%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera

Separation 5% (deselected), Parallax Shift 0% (deselected)

3D Image Processing - Panel B

3D Volume Projection = Precise

Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels

Background = Black

Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 35%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera

Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Panel A shows a three-dimensional volume projection of PARP1 Alexa Fluor™ 532 Direct Conjugate with DAPI, acquired by lattice structured illumination microscopy (SIM). Panel B presents a single 2D optical section (Z-plane 56 of 142); the white box denotes the region of interest used to generate the corresponding volume projection. PARP1 exhibits robust nuclear localization consistent with prior confocal and Airyscan analyses, with markedly increased spatial definition at lattice SIM resolution (~60-80 nm in the X-Y plane). At this resolution, PARP1 aggregates associated with DNA damage-induced PARylation become more clearly resolved, though they remain partially unresolved relative to single-molecule modalities. Notably, PARylation is a highly dynamic process; the observed aggregate morphologies likely represent different functional states dependent on damage context and timing, and definitive assignment of PARP1 structural states would require a dedicated temporal analysis.

Image Corrections

NONE

Image by P. Raut


Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000053
Stage	SIM
Instrument	ZEISS Elyra
Microscope configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Structured source metadata values	65
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	BA002CC3AD2CF4D9A09D7E69F79A55ADEC4110EEEEAFB15A39AB84460CBCD28A
Method PDF SHA-256	78B7227238B9019A8D4461CFCF41645631EE66112AEBC3AAEC42353B32B03EBF

ORIGINAL IMAGE EVIDENCE / SIM²



EVIDENCE TIER	SIM ²
INSTRUMENT	ZEISS Elyra
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 532; 647
METADATA VALUES	67

STRUCTURED ACQUISITION METADATA / SIM²

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	110 Slices (3.27 µm)
Channels	3
Scaling (Per Pixel)	0.02111 µm X 0.02111 µm X 0.03000 µm
Image size (Pixels)	441 X 358 192 X 218
Image Size (Scaled)	9.31 µm X 7.56 µm 4.05 µm X 4.60 µm
Bit Depth	16 Bit
Image Center Position	X: 62.69 µm, Y: 35.57 µm
Stage Position	X: 62.69 µm, Y: 35.57 µm
ROI Center Offset	X: 0.00 µm, Y: 0.00 µm
Reflectors	SR HE LBF 405/488/561/642
Beam Splitter	405, 488, 561, 642
Filter Excitation Wavelength	200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method	Fluorescence
Illumination Wavelength	488 640 405
Channel Name	T2-Axiocam 820 #2-SR T1-Axiocam 820 -SR T3-Axiocam 820 #2-SR
Excitation Wavelength	488 640 405
Emission Wavelength	525 729 445
Effective NA	1.4
Detection Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800 503 - 546, 657 - 800
Imaging Device	Axiocam 820 #2 Axiocam 820
Camera Adapter	1X
Exposure Time	150 ms 100 ms 50 ms
Depth of Focus	0.81 µm 1.13 µm 0.69 µm
Binning Mode	2,2
Laser Wavelength	488 nm @4.0% 640 nm @2.0% 405 nm @2.0%
Frame Time	1.59 s 1.33 s
Scan Direction	Unidirectional
Grating	36.5 µm grating 27.5 µm grating
Result Sampling	4
Process Sampling	4

FIELD	SOURCE-DOCUMENTED VALUE
Algorithm	SIM2
Input SNR	Low
Regularization Weight	0.0000
Iterations	3
Filter	No Filter
Fitting	Fast Fit
Normalization	Auto
3D Volume Projection	Precise
Appearance	Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels Threshold 20%, Ambient Light 0%, Specular Light 100%, Shininess 35%
Background	Black
Light	Brightness 125%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected) View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
3D Surface Projection	Precise

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - (BLR018M)
Cat ID# - DC 000053

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L)
Cat ID# - SA 000059

PARP1 (Poly(ADP-ribose) polymerase 1) has a multifaceted function in the detection of DNA damage, including single-strand breaks (SSBs) and double-strand breaks (DSBs) via binding at the site of damage. PARP1 catalyzes the synthesis of poly(ADP-ribose) (PAR) chains using NAD⁺; the resulting PAR chains signal the recruitment of DNA repair proteins to the break, initiating the repair process. PARP1, as well as PARP2, are retained in nucleoli via interaction with the multifunctional nucleolar hub protein nucleophosmin, which is implicated in multiple steps of ribosome biogenesis, including rDNA transcription and elongation, as well as rRNA processing.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Rabbit anti-Human PARP1 Recombinant Monoclonal - Identifyn® SRM Alexa Fluor™ 532 - Direct Conjugate, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS.

The second primary antibody, Nuclear Pore Complex Monoclonal Antibody (39C7), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Nuclear Pore Complex(NPC), Identifyn® SRM Alexa Fluor™ 647 Staining not represented in data analysis

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:

Z Stack - (3D) - Panel A

Z-Stack = 110 Slices (3.27 µm)

Channels = 3
Scaling (Per Pixel) = 21.11 nm X 21.11 nm X 30.00 nm
Image size (Pixels) = 441 X 358
Image Size (Scaled) = 9.31 µm X 7.56 µm
Bit Depth = 16 Bit
Image Center Position = X: 62.69 µm, Y: 35.57 µm
Stage Position = X: 62.69 µm, Y: 35.57 µm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D) - Panel B

Z-Stack = 110 Slices (3.27 µm)

Channels = 3
Scaling (Per Pixel) = 21.11 nm X 21.11 nm X 30.00 nm
Image size (Pixels) = 192 X 218
Image Size (Scaled) = 4.05 µm X 4.60 µm
Bit Depth = 16 Bit
Image Center Position = X: 62.69 µm, Y: 35.57 µm
Stage Position = X: 62.69 µm, Y: 35.57 µm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

532 -

Reflectors = SR HE LBF 405/488/561/642
Beam Splitter = 405, 488, 561, 642
Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method = Fluorescence
Illumination Wavelength = 488
Channel Name = T2-Axiocam 820 #2-SR
Excitation Wavelength = 488
Emission Wavelength = 525
Effective NA = 1.4
Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Imaging Device = Axiocam 820 #2
Camera Adapter = 1X
Exposure Time = 150 ms
Depth of Focus = 0.81 µm
Binning Mode = 2,2
Laser Wavelength = 488 nm @4.0%
Frame Time = 1.59 s
Scan Direction = Unidirectional
Grating = 36.5 µm grating

647 - (Not Shown in Data Analysis)

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 640
 Channel Name = T1-Axiocam 820 -SR
 Excitation Wavelength = 640
 Emission Wavelength = 729
 Effective NA = 1.4
 Detection Wavelength = 503 - 546, 657 - 800
 Imaging Device = Axiocam 820
 Camera Adapter = 1X
 Exposure Time = 100 ms
 Depth of Focus = 1.13 μm
 Binning Mode = 2,2
 Laser Wavelength = 640 nm @2.0%
 Frame Time = 1.33 s
 Scan Direction = Unidirectional
 Grating = 27.5 μm grating

DAPI -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 405
 Channel Name = T3-Axiocam 820 #2-SR
 Excitation Wavelength = 405
 Emission Wavelength = 445
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820 #2
 Camera Adapter = 1X
 Exposure Time = 50 ms
 Depth of Focus = 0.69 μm
 Binning Mode = 2,2
 Laser Wavelength = 405 nm @2.0%
 Frame Time = 1.33 s
 Scan Direction = Unidirectional
 Grating = 36.5 μm grating

Post Processing - SIM2 Processing

Processing Dimension - 3D
 Result Sampling = 4
 Process Sampling = 4
 Channels = 3
 Algorithm = SIM2
 Input SNR = Low
 Regularization Weight = 0.0000
 Iterations = 3
 Filter = No Filter
 Fitting = Fast Fit
 Normalization = Auto

3D Image Processing - Panel A

3D Volume Projection = Precise
 Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 125%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

3D Image Processing - Panel B

3D Surface Projection = Precise
 Appearance = Threshold 20%, Ambient Light 0%, Specular Light 100%, Shininess 35%
 Background = Black
 Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Panel A shows a three-dimensional volume projection of PARP1 detected with the Alexa Fluor™ 532 direct conjugate (white) and DAPI (blue). A white box denotes the region of interest (ROI) analyzed in Panel B, which presents a three-dimensional surface projection of the same PARP1 assemblies. At SIM² lateral resolution (~40-60 nm), PARP1 displays nuclear localization and aggregate morphologies comparable to those observed by lattice SIM, consistent with the relatively large size of PARP1 assemblies. Structural differentiation between SIM and SIM² is therefore limited at this stage; further resolution of PARP1 sub-units and higher-order organization requires single-molecule localization microscopy.

Image Corrections

NONE

Image by P. Raut

Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000053
Stage	SIM²
Instrument	ZEISS Elyra
Microscope configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Structured source metadata values	67
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	3BE4A56260F5238D532F0C6516DCFDE10EAB3A24CA57CC1FE9A9E599D513995F
Method PDF SHA-256	E1FA0A265E01EC4B61CF948747034F1139D7303B4EF56F7C4D0E6C0846AE1419