

IDENTIFYN RESOLUTION STEPDOWN

H2AX

MICROSCOPY EVIDENCE ASSESSMENT

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647

WESTERN BLOT -> WIDEFIELD -> WIDEFIELD SIM / APOTOME -> CONFOCAL -> AIRYSCAN
SIM -> SIM² -> SINGLE-MOLECULE STORM

8


ORDERED EVIDENCE TIERS

0

SUPPORTING CONTROL

PENDING

SCIENTIST REVIEW



REPRESENTATIVE HIGHEST-AVAILABLE VIEW / SINGLE-MOLECULE STORM

DC-000035 | PRODUCT-LEVEL EVIDENCE DOSSIER

SOURCE-GROUNDED MICROSCOPY + BIOLOGICAL REVIEW / IDENTIFYN EVIDENCE PRESERVED

CRO CAPABILITY DEMONSTRATED

Identifyn performed the experimental and imaging work documented below. BENNETT converts that operating history into a governed, source-bounded scientific record.

- Source-documented sample preparation and antibody/reagent workflow
- Western blot imaging
- Widefield fluorescence microscopy
- Widefield structured illumination / Apotome optical sectioning
- Confocal microscopy
- Airyscan acquisition and reconstruction workflow
- Structured illumination microscopy
- SIM² reconstruction workflow
- Single-molecule STORM presentation workflow
- Treatment-response experimental workflow
- Z-stack / 3D acquisition
- Multi-channel imaging

SCIENTIFIC QUESTION

What source-bounded experimental and microscopy evidence is preserved for H2AX, and what is the strongest interpretation allowed by the available modalities and controls?

EVIDENCE AVAILABLE

Western blot -> Widefield -> Widefield SIM — Apotome -> Confocal -> Airyscan -> SIM -> SIM² -> Single-molecule STORM

BENNETT ASSESSMENT

The preserved DC-000035 record contains the full ordered evidence sequence from Western blot through single-molecule STORM. Each modality retains its own acquisition and processing boundary.

BENNETT / MICROSCOPY AND BIOLOGICAL EVIDENCE REVIEW

SOURCE-GROUNDED RETROSPECTIVE REVIEW / PENDING SCIENTIST REVIEW

SOURCE STATE	SOURCE HASH VERIFIED / EVENT RECORDED
BENNETT EVIDENCE REVIEW	COMPLETED / SOURCE-BOUNDED
NATIVE IMAGE REPROCESSING	NOT PERFORMED
LITERATURE SOURCES	INDEPENDENTLY VERIFIED
SCIENTIST REVIEW	PENDING
PUBLICATION	NOT AUTHORIZED

OVERALL CALL / FULL STEPDOWN / SOURCE-BOUNDED REVIEW

The preserved DC-000035 record contains the full ordered evidence sequence from Western blot through single-molecule STORM. Each modality retains its own acquisition and processing boundary.

SETTINGS ASSESSMENT

The record preserves 8 ordered evidence tier(s). Structured source metadata values are reported by stage and remain tied to the original method objects.

BIOLOGICAL SUPPORT

The record supports retrospective, source-bounded microscopy review. Target identity, specificity, treatment effect, binding, interaction and mechanism are not inferred from presentation images alone.

CONTROL ASSESSMENT

No separate control record is preserved; none is inferred from method or template language.

CLAIM CEILING

This draft supports historical capability documentation, source registration, modality ordering, settings review and bounded interpretation. It does not independently validate antibody specificity, quantify a population response, establish binding, interaction, causation or mechanism, or authorize publication.

PROCESSING DISCLOSURE

Retrospective BENNETT architectural evaluation of preserved Identifyn source evidence. BENNETT did not perform native-file reprocessing, new deterministic measurement or the historical image post-processing represented in this record.

REAGENTS AND ACQUIRED CHANNELS

Preparation reagents are reported separately from channels actually documented in the acquisition settings. Fluorophores mentioned in staining are not automatically treated as acquired channels.

EVIDENCE OBJECT	REAGENTS / FLUOROPHORES PRESENT	CHANNELS ACTUALLY ACQUIRED
Western blot	647	Not deterministically stated in the acquisition block
Widefield	405/DAPI; 532; 555; 647	3; 50 Cy 5; AF532-49014; 49 DAPI; 625 - 655; 515 - 545; 335 - 383; 665 - 715
Widefield SIM — Apotome	405/DAPI; 532; 555; 647	3; 50 Cy 5; AF532-49014; 49 DAPI; 625 - 655; 515 - 545; 335 - 383; 665 - 715
Confocal	405/DAPI; 488; 532; 647	3; None; 639 nm @0.2%; 488 nm @0.05%; 405 nm @0.2%; 653; 528; 353
Airyscan	405/DAPI; 488; 532; 647; 750	3; None; 639 nm @0.4%; 488 nm @0.5%; 405 nm @0.2%; 653; 528; 353
SIM	405/DAPI; 488; 647	3; SR HE LBF 405/488/561/642; 200 - 405, 200 - 488, 200 - 561, 200 - 642; 419 - 470, 503 - 546, 576 - 617, 657 - 800; T1-Axiocam 820-SR; T2-Axiocam 820-SR; T3-Axiocam 820 #2-SR; 640
SIM ²	405/DAPI; 488; 647	3; SR HE LBF 405/488/561/642; 200 - 405, 200 - 488, 200 - 561, 200 - 642; 419 - 470, 503 - 546, 576 - 617, 657 - 800; T1-Axiocam 820-SR; T2-Axiocam 820-SR; T3-Axiocam 820 #2-SR; 640
Single-molecule STORM	647	Both Channels

MODALITY ASSESSMENT / PART 1 OF 8

Western blot / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: Azure Sapphire FL. Objective: not explicitly stated. Platform: Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s). Structured source metadata values: 11.

VISIBLE OBSERVATION

A preserved presentation image is available for Western blot. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Western blot, documents platform context as Azure Sapphire FL, and records detected dye/channel descriptors as 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

Band position and displayed intensity do not independently establish analytical specificity, target identity, linearity, or treatment effect.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Western blot evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 2 OF 8

Widefield / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Axio Imager.M1. Objective: Plan-Apochromat 63x/1.40 Oil M27. Platform: Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705, and X-Cite Xylis Lamp, was used with the following parameters for each dye: Acquisition: Channels: 3; Scaling (Per Pixel): 0.110 μm X 0.110 μm ; Image size (Pixels): 1216 X 1028; Image Size (Pixels): 1216 X 1028; Bit Depth: 14 Bit. Detector/camera: Imaging Device: Axiocam 705. Timing: Exposure Time: 200 ms; 100 ms. Illumination: Light source: X-Cite Xylis Lamp Intensity @12.00%; X-Cite Xylis Lamp Intensity @10.00%. Optical/scan: Effective NA: 1.4. Structured source metadata values: 40.

VISIBLE OBSERVATION

A preserved presentation image is available for Widefield. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Widefield, documents platform context as ZEISS Axio Imager.M1, and records detected dye/channel descriptors as 405/DAPI; 532; 555; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

A presentation image supports gross localization and source-bounded co-occurrence, not direct binding, molecular colocalization, or native quantitative measurement.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Widefield evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 3 OF 8

Widefield SIM — Apotome / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Axio Observer.Z1/7 Apotome. Objective: EC Plan-Neofluar 63X/1.25 oil M27. Platform: Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:. Acquisition: Z-Stack: 86 Slices (8.5 μm); Channels: 3; Scaling (Per Pixel): 0.143 μm X 0.143 μm X 0.100 μm ; Image size (Pixels): 786 X 689; Image Size (Pixels): 786 X 689; Bit Depth: 14 Bit. Detector/camera: Imaging Device: 807 mono. Timing: Exposure Time: 200 ms; 50 ms. Illumination: Light source: X-Cite Xylis Lamp Intensity @15.00%; X-Cite Xylis Lamp Intensity @8.00%. Optical/scan: Effective NA: 1.25. Processing: Projection: View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 47.

VISIBLE OBSERVATION

A preserved presentation image is available for Widefield SIM — Apotome. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Widefield SIM — Apotome, documents platform context as ZEISS Axio Observer.Z1/7 Apotome, and records detected dye/channel descriptors as 405/DAPI; 532; 555; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

Improved optical sectioning does not by itself establish reconstructed super-resolution or molecular interaction.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Widefield SIM — Apotome evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 4 OF 8

Confocal / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS LSM 980. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:. Acquisition: Z-Stack: 181 Slices (5.4 μm); Channels: 3; Scaling (Per Pixel): 0.059 μm X 0.059 μm X 0.030 μm ; Image size (Pixels): 1092 X 728; Image Size (Pixels): 1092 X 728; Bit Depth: 16 Bit. Detector/camera: Imaging Device: GaAsP-Pmt3; Detector Type: GaAsP-PMT. Timing: Pixel Time: 0.30 μs ; Frame Time: 3.69 s. Illumination: Laser Wavelength: 639 nm @0.2%; 488 nm @0.05%. Optical/scan: Pinhole: 1.00 AU / 66 μm ; 1.00 AU / 53 μm ; Scan Zoom: 2.0; LSM Scan Speed: 12; Scan Direction: Bidirectional; Averaging: 8; Effective NA: 1.4. Processing: Projection: View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black; LSM Plus Mode: 7.2 (3D, Auto); 7.8 (3D, Auto). Structured source metadata values: 55.

VISIBLE OBSERVATION

A preserved presentation image is available for Confocal. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Confocal, documents platform context as ZEISS LSM 980, and records detected dye/channel descriptors as 405/DAPI; 488; 532; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy and antibody-validation guidance that presentation evidence must remain bounded by modality, controls and source documentation.

INTERPRETIVE LIMIT

Optical sectioning and source-bounded 3D presentation do not establish binding, mechanism, or population-level effect.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Confocal evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 5 OF 8

Airyscan / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS LSM 980. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:. Acquisition: Z-Stack: 181 Slices (5.4 μm); Channels: 3; Scaling (Per Pixel): 35.29 nm X 35.29 nm X 30.00 nm; Image size (Pixels): 1387 X 1105; 679 X 653; Image Size (Pixels): 1387 X 1105; 679 X 653; Bit Depth: 16 Bit. Detector/camera: Imaging Device: Airyscan; Detector Type: GaAsP-PMT. Timing: Pixel Time: 1.21 μs ; Frame Time: 17.06 s. Illumination: Laser Wavelength: 639 nm @0.4%; 488 nm @0.5%. Optical/scan: Pinhole: 5.00 AU / 328 μm ; 5.00 AU / 265 μm ; Scan Zoom: 2.2; LSM Scan Speed: 6; Scan Direction: Bidirectional; Averaging: 4; Effective NA: 1.4. Processing: Projection: View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black; AiryScan Mode: SuperResolution 9.6 (3D, Auto); SuperResolution 10.0 (3D, Auto). Structured source metadata values: 63.

VISIBLE OBSERVATION

A preserved presentation image is available for Airyscan. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Airyscan, documents platform context as ZEISS LSM 980, and records detected dye/channel descriptors as 405/DAPI; 488; 532; 647; 750.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Displayed feature separation remains reconstruction dependent; no unsupported numerical resolution or molecular architecture claim is permitted.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Airyscan evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 6 OF 8

SIM / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Elyra. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye: Acquisition: Z-Stack: 124 Slices (3.69 μ m); Channels: 3; Scaling (Per Pixel): 16.89 nm X 16.89 nm X 30.00 nm; Image size (Pixels): 5056 X 5056; 2010 X 2426; Image Size (Pixels): 5056 X 5056; 2010 X 2426; Bit Depth: 16 Bit. Detector/camera: Imaging Device: AxioCam 820; AxioCam 820 #2-SR. Timing: Exposure Time: 120 ms; 100 ms; Frame Time: 1.59 s; 1.33 s. Illumination: Laser Wavelength: 640 nm @3.00%; 488 nm @0.9%; Illumination Wavelength: 640; 488. Optical/scan: Scan Direction: Unidirectional; Effective NA: 1.4. Processing: Algorithm: SIM; Process Sampling: N/A; Result Sampling: 4; Projection: View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 70.

VISIBLE OBSERVATION

A preserved presentation image is available for SIM. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as SIM, documents platform context as ZEISS Elyra, and records detected dye/channel descriptors as 405/DAPI; 488; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Reconstruction-dependent structural separation requires native reconstruction QC before numerical resolution or molecular dimensions can be claimed.

CLAIM CEILING

This tier supports a source-bounded statement that the historical SIM evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 7 OF 8

SIM² / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: ZEISS Elyra. Objective: Plan-Apochromat 63X/1.40 oil DIC M27. Platform: Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye: Acquisition: Z-Stack: 154 Slices (4.59 μ m); Channels: 3; Scaling (Per Pixel): 16.89 nm X 16.89 nm X 30.00 nm; Image size (Pixels): 5056 X 5056; 2290 X 1950; Image Size (Pixels): 5056 X 5056; 2290 X 1950; Bit Depth: 16 Bit. Detector/camera: Imaging Device: AxioCam 820; AxioCam 820 #2. Timing: Exposure Time: 120 ms; 100 ms; Frame Time: 1.59 s; 1.33 s. Illumination: Laser Wavelength: 640 nm @3.00%; 488 nm @0.9%; Illumination Wavelength: 640; 488. Optical/scan: Scan Direction: Unidirectional; Effective NA: 1.4. Processing: Algorithm: SIM²; Process Sampling: 4; Result Sampling: 4; Projection: View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); View Angle 19°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected); Background: Black. Structured source metadata values: 71.

VISIBLE OBSERVATION

A preserved presentation image is available for SIM². BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as SIM², documents platform context as ZEISS Elyra, and records detected dye/channel descriptors as 405/DAPI; 488; 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Further reconstruction-dependent segmentation does not establish molecular dimensions, structure identity, or mechanism without independent support.

CLAIM CEILING

This tier supports a source-bounded statement that the historical SIM² evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

MODALITY ASSESSMENT / PART 8 OF 8

Single-molecule STORM / SOURCE EVIDENCE PRESERVED

DOCUMENTED SETTINGS

Platform: Bruker Vutara VXL. Objective: Olympus 60X/1.30 silicon oil objective. Platform: Bruker Vutara VXL Biplane Single Molecule Microscope, using an Olympus 60X/1.30 silicon oil objective, was used with the following parameters:.

Detector/camera: Camera: Both Cameras. Timing: Exposure Time: 20ms Selected; Super Resolution Camera Exposure Time: 20ms; Widefield Camera Exposure Time: 100ms. Illumination: Photoactivation/Imaging Laser: 0% Selected.

Structured source metadata values: 80.

VISIBLE OBSERVATION

A preserved presentation image is available for Single-molecule STORM. BENNETT records the visible source presentation without creating a new native measurement.

SOURCE-REPORTED RESULT

The source identifies this object as Single-molecule STORM, documents platform context as Bruker Vutara VXL, and records detected dye/channel descriptors as 647.

LITERATURE CORRELATION

Consistent with peer-reviewed microscopy guidance that advanced-resolution outputs remain dependent on acquisition, reconstruction or localization quality control.

INTERPRETIVE LIMIT

Without localization tables, uncertainty, blinking, drift, and independent cluster evidence, the presentation cannot establish molecule count, cluster truth, exact precision, or molecular architecture.

CLAIM CEILING

This tier supports a source-bounded statement that the historical Single-molecule STORM evidence and documented settings are preserved for review; stronger biological or quantitative conclusions require compatible native data, controls and scientist review.

CROSS-MODALITY SYNTHESIS

The preserved DC-000035 record contains the full ordered evidence sequence from Western blot through single-molecule STORM. Each modality retains its own acquisition and processing boundary.

Different acquisitions are not treated as a quantitative intensity ladder. Any continuity statement remains qualitative and source bounded.

BENNETT LITERATURE CORRELATION — SOURCES INDEPENDENTLY VERIFIED

BENNETT uses independently verified primary scientific sources to test whether the interpretation is scientifically consistent and properly bounded. Literature correlation does not retroactively validate the Identifyn dataset or replace scientist review.

- Antibody and microscopy evidence require orthogonal controls and modality-appropriate interpretation; presentation images alone do not establish analytical specificity.
- Advanced-resolution presentation depends on acquisition, reconstruction or localization quality control before numerical resolution, molecular dimensions or cluster identity can be claimed.

REFERENCES

1. Uhlen M, et al. A proposal for validation of antibodies. Nat Methods. 2016;13:823-827. doi:10.1038/nmeth.3995.
2. Gustafsson MGL. Surpassing the lateral resolution limit by a factor of two using structured illumination microscopy. J Microsc. 2000;198:82-87. doi:10.1046/j.1365-2818.2000.00710.x.
3. Sage D, et al. Quantitative evaluation of software packages for single-molecule localization microscopy. Nat Methods. 2015;12:717-724. doi:10.1038/nmeth.3442.

RELEASE BOUNDARY

This assessment remains a draft pending scientist review. No publication is authorized. Any future native-file analysis must enter BENNETT as a new evidence snapshot with routed engine authority, native-data QA, methods, limitations, and an independently reviewed claim ceiling.

SOURCE RECONCILIATION / QA / PROVENANCE

SOURCE STATUS - PRESERVED / AUTO-RECONCILED

Historical Identifyn source material was subjected to BENNETT retrospective QA. Original source state and correction history remain preserved in provenance.

HIGH CONFIDENCE - AUTO-RECONCILED

Airyscan / Scaling (Per Pixel): 35.29 nm X 35.29 nm X 30.00 nm -> 0.03529 μ m X 0.03529 μ m X 0.03000 μ m. Both numerical magnitude and unit were normalized, then checked against the reported field dimensions; no unit label was changed alone. Supporting evidence: Image Size (Pixels)=1387 X 1105; Image Size (Scaled)=48.94 μ m X 38.99 μ m; PASS - INTERNALLY CONSISTENT. Documented sampling reproduces the reported physical field dimensions. Maximum field-dimension error 0.01%.

HIGH CONFIDENCE - AUTO-RECONCILED

SIM / Scaling (Per Pixel): 16.89 nm X 16.89 nm X 30.00 nm -> 0.01689 μ m X 0.01689 μ m X 0.03000 μ m. Both numerical magnitude and unit were normalized, then checked against the reported field dimensions; no unit label was changed alone. Supporting evidence: Image Size (Pixels)=5056 X 5056; Image Size (Scaled)=85.40 μ m X 85.40 μ m; PASS - INTERNALLY CONSISTENT. Documented sampling reproduces the reported physical field dimensions. Maximum field-dimension error 0.00%.

HIGH CONFIDENCE - AUTO-RECONCILED

SIM² / Scaling (Per Pixel): 16.89 nm X 16.89 nm X 30.00 nm -> 0.01689 μ m X 0.01689 μ m X 0.03000 μ m. Both numerical magnitude and unit were normalized, then checked against the reported field dimensions; no unit label was changed alone. Supporting evidence: Image Size (Pixels)=5056 X 5056; Image Size (Scaled)=85.40 μ m X 85.40 μ m; PASS - INTERNALLY CONSISTENT. Documented sampling reproduces the reported physical field dimensions. Maximum field-dimension error 0.00%.

HIGH CONFIDENCE - SCIENTIST-AUTHORIZED CORRECTION

Cross-tier / treatment concentration: 200 μ M; 2 μ M -> Etoposide = 2 μ M. The governing scientist identified every 200 μ M occurrence as a technician transcription error. The canonical historical concentration is 2 μ M. Supporting evidence: Scientist-authorized correction supplied for the Identifyn historical archive; original technician text remains preserved unchanged. Authority: Dr. Brian T. Bennett. Preservation: Original source text and hashes remain unchanged; only the derivative canonical field is corrected.

UNRESOLVED - PRESERVED AS CONFLICT

Cross-tier / treatment duration: {"Western blot": ["24 hours"], "Single-molecule STORM": ["6 hours"]} -> Preserved as modality-specific historical duration values. Genuinely different treatment durations are present across evidence stages, and no source or scientist correction resolves them. Supporting evidence: Only duration values from explicit etoposide-treatment sentences were classified; concentration, wavelength, resolution and pixel-size values were excluded.

BENNETT EVIDENCE REVIEW COMPLETE

INTERPRETATION TO EVIDENCE

BENNETT has completed the bounded retrospective assessment.
The following pages switch ownership to the preserved Identifyn record.



IDENTIFYN PRESERVED SOURCE EVIDENCE

ORIGINAL RECORD BEGINS

Images, source methods, acquisition settings and provenance follow in the locked
Stepdown order. Scientific wording is preserved; missing modalities are not invented.
No separate control record is present in the preserved product record.

PRODUCT-LEVEL STEPDOWN MAP


This dossier preserves the complete available evidence progression for DC-000035. Each tier remains tied to its original image, source method, instrument context and cryptographic identity. No separate control record is present in the preserved product record, and none is inferred.

CANONICAL STEPDOWN	EVIDENCE STAGE	INSTRUMENT	STATUS
1	Western blot	Azure Sapphire FL	PRESERVED
2	Widefield	ZEISS Axio Imager.M1	PRESERVED
3	Widefield SIM - Apotome	ZEISS Axio Observer.Z1/7 Apotome	PRESERVED
4	Confocal	ZEISS LSM 980	PRESERVED
5	Airyscan	ZEISS LSM 980	PRESERVED
6	SIM	ZEISS Elyra	PRESERVED
7	SIM ²	ZEISS Elyra	PRESERVED
8	Single-molecule STORM	Bruker Vutara VXL	PRESERVED

INTERPRETIVE BOUNDARY

Ordering evidence by resolution tier does not make the modalities interchangeable. Each stage retains its acquisition environment and scientific limitations. This pilot organizes the record; it does not synthesize new biological claims.

ORIGINAL IMAGE EVIDENCE / WESTERN BLOT



EVIDENCE TIER	Western blot
INSTRUMENT	Azure Sapphire FL
CELL MODEL	HeLa
DETECTED DYES	647
METADATA VALUES	11

STRUCTURED ACQUISITION METADATA / WESTERN BLOT

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)
Scanner	Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)
Detection mode	Fluorescence
Laser	638nm
Emission Filter	676±37nm
Detector	PMT
Intensity	2
Pixel size	100 µm
Scan Speed	High
Quality	High
Focus	Membrane (0 µm in Z-axis)
Image	K. Small

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Expected Molecular Weight: 15kDa

HeLa cells were damaged with 200 µM etoposide for 24 hours. Cells were then homogenized using an ultrasonication probe and RIPA lysis buffer, and then centrifuged to pellet the cell debris. The lysate supernatant (100 µg) was mixed with NuPAGE™ LDS Sample Buffer (4X), heated at 80 °C for 10 min, and subjected to electrophoresis using NuPAGE™ 4 to 12% Bis-Tris, 1.0-1.5 mm, Mini Protein Gels in a Mini Gel Tank. The ThermoFisher's PageRuler™ Prestained Protein Ladder, 10 to 180 kDa, was used as a molecular weight marker.

Post electrophoresis, the proteins were transferred to a polyvinylidene difluoride (PVDF) membrane using the iBlot™ 3 Starter Kit, PVDF. The transfer membrane was blocked with Blocker™ Casein for 1 hour and then washed five times with PBS. The blocked membrane was incubated with Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647 - Direct Conjugate for 2 hours, followed by five washes in 1X PBS. The membrane was dried for 10 minutes at room temperature before imaging.

Scanner = Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)

Scan Settings

Detection mode = Fluorescence

Laser = 638nm

Emission Filter = 676±37nm

Detector = PMT

Intensity = 2

Pixel size = 100 µm

Scan Speed = High

Quality = High

Focus = Membrane (0 µm in Z-axis)

Post Processing

None


Image: K. Small

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	Western blot
Instrument	Azure Sapphire FL
Microscope configuration	Azure Sapphire FL with Customizable Laser and Filter Modules, Enabling UV- NIR and Z-Plane(s)
Structured source metadata values	11
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	862B2895D2D6697CA80EB15B603CA694C958670EE1ECA5F2E7336F2927056793
Method PDF SHA-256	4588D70073037AB76ADB882FC73B7D3EDFC324C02406977F97EB28543ADE3F2A

ORIGINAL IMAGE EVIDENCE / WIDEFIELD



EVIDENCE TIER	Widefield
INSTRUMENT	ZEISS Axio Imager.M1
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 532; 555; 647
METADATA VALUES	40

STRUCTURED ACQUISITION METADATA / WIDEFIELD

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705, and X-Cite Xylis Lamp, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63x/1.40 Oil M27
Channels	3
Scaling (Per Pixel)	0.110 μ m X 0.110 μ m
Image size (Pixels)	1216 X 1028
Image Size (Scaled)	133.18 μ m X 112.59 μ m
Bit Depth	14 Bit
Image Center Position	X: 0.00 μ m, Y: 0.00 μ m
ROI Center Offset	X: 0.00 μ m, Y: 0.00 μ m
Reflector	50 Cy 5 AF532-49014 49 DAPI
Beam Splitter	660 550 395
Filter Excitation Wavelength	625 - 655 515 - 545 335 - 383
Filter Emission Wavelength	665 - 715 555 - 595 420 - 470
Contrast Method	Fluorescence
Light source	X-Cite Xylis Lamp Intensity @12.00% X-Cite Xylis Lamp Intensity @10.00% X-Cite Xylis Lamp Intensity @15.00%
Exposure Time	200 ms 100 ms 20 ms
Excitation Wavelength	653 528 353
Emission Wavelength	668 553 465
Effective NA	1.4
Imaging Device	Axiocam 705
Camera Adapter	1X
Depth of Focus	1.03 μ m 0.86 μ m 0.72 μ m
Binning Mode	2,2
Profile Definition	Stroke Thickness 1, Normal View, Grid Distance 1 Profile View = X and Y axes were set to Auto - X (Min 0.0 and Max 90.3), Y (Min 28 and Max 1256)

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Identifyn® Secondary Antibody

Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat Anti-Rabbit, IgG (H + L)

Cat ID# - SA 000024

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ -H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ -H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ -H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ -H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with Thermo Fisher, Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS.

The second primary antibody, Golgi GOLPH2 Polyclonal Antibody, was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat

Anti-Rabbit, IgG (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI, cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, Axiocam 705, and X-Cite Xylis Lamp, was used with the following parameters for each dye:

Single Plane (2D)

Channels = 3

Scaling (Per Pixel) = 0.110 μm X 0.110 μm

Image size (Pixels) = 1216 X 1028

Image Size (Scaled) = 133.18 μm X 112.59 μm

Bit Depth = 14 Bit

Image Center Position = X: 0.00 μm , Y: 0.00 μm

ROI Center Offset = X: 0.00 μm , Y: 0.00 μm

647 -

Reflector = 50 Cy 5

Beam Splitter = 660

Filter Excitation Wavelength = 625 - 655

Filter Emission Wavelength = 665 - 715

Contrast Method = Fluorescence

Light source = X-Cite Xylis Lamp Intensity @12.00%

Exposure Time = 200 ms

Excitation Wavelength = 653

Emission Wavelength = 668

Effective NA = 1.4

Imaging Device = Axiocam 705

Camera Adapter = 1X

Depth of Focus = 1.03 μm

Binning Mode = 2,2

532 -

Reflector = AF532-49014
 Beam Splitter = 550
 Filter Excitation Wavelength = 515 - 545
 Filter Emission Wavelength = 555 - 595
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @10.00%
 Exposure Time = 100 ms
 Excitation Wavelength = 528
 Emission Wavelength = 553
 Effective NA = 1.4
 Imaging Device = Axiocam 705
 Camera Adapter = 1X
 Depth of Focus = 0.86 μ m
 Binning Mode = 2,2

DAPI -

Reflector = 49 DAPI
 Beam Splitter = 395
 Filter Excitation Wavelength = 335 - 383
 Filter Emission Wavelength = 420 - 470
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @15.00%
 Exposure Time = 20 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Imaging Device = Axiocam 705
 Camera Adapter = 1X
 Depth of Focus = 0.72 μ m
 Binning Mode = 2,2

Split View - Panel A

The top-left panel features Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), the top-right panel shows Golgi visualized with Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat Anti-Rabbit, IgG (H + L), the bottom-left panel features DAPI nuclear staining, and the bottom-right panel features the merged image of all channels.

Profile View Display - Panel B & C

Profile Definition = Stroke Thickness 1, Normal View, Grid Distance 1 Profile View = X and Y axes were set to Auto - X (Min 0.0 and Max 90.3), Y (Min 28 and Max 1256)

Summary

Panel A- A split view illustrates the signal and localization of Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647 IgG1, direct conjugate. In the upper left, the image shows the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647 IgG1, direct conjugate (red). The upper right displays the Golgi primary antibody labeled with Identifyn® SRM Alexa Fluor™ 532 (yellow). The lower left features DAPI nuclear staining (blue). Finally,

the lower right presents a merged image combining all three channels.

Note - We observe that the localization of the H2AX/647 direct conjugate is nuclear with a "pan" distribution, similar to what is seen with the unconjugated primary antibody. This pan-nuclear staining is expected at diffraction-limited widefield resolution; however, as we shift towards super-resolution imaging, the sample will begin to reveal structural details.

Panels B and C - These panels present a profile view, where the image in Panel C is next to the histogram in Panel B. A line has been drawn through several cells, specifically through the nuclei. This line measures both the distance and fluorescent intensity. We use this measurement to display the intensity values of the Mouse anti-Human H2AxMonoclonal - Identifyn® SRM Alexa Fluor™ 647 IgG1, direct conjugate (red) at an LED excitation power of 12% and an exposure time of 200ms, thus confirming the localization and high brightness of the direct conjugate. Additionally, we include the Golgi primary antibody, counterstained with Identifyn® SRM Alexa Fluor™ 532, in the measurement.

Image Correction

None

Image by E.F. Simon

Image Analysis by B.T. Bennett


Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	Widefield
Instrument	ZEISS Axio Imager.M1
Microscope configuration	Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.40 Oil M27, AxioCam 705, and X-Cite Xylis Lamp, was used with the following parameters for each dye:
Structured source metadata values	40
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	510EDE1EBA63B9BBC92E8D1E632507513F8A31C0D01BD258C352A884921FF917
Method PDF SHA-256	F6D106A9C41823A0BD045EED6A6C7BDAD49B169A7CE6C8996F6A82D85273BDAB

ORIGINAL IMAGE EVIDENCE / WIDEFIELD SIM - APOTOME



Panel A – This panel presents a three-dimensional volume projection of the Mouse anti-Human H2Ax Monoclonal - Identifyn™ SRM Alexa Fluor™ 647 (red) conjugate and a Golgi primary antibody labeled with Identifyn™ SRM Alexa Fluor™ 532 (yellow), utilizing widefield structured illumination microscopy (SIM) resolution.

Panel B – This panel displays a graphical colocalization view featuring the Mouse anti-Human H2Ax Monoclonal - Identifyn™ SRM Alexa Fluor™ 647 (X-axis) and DAPI (Y-axis).

Panel C – This panel shows the image from which the colocalization measurements were obtained at z-plane 47 out of a total of 86 collected planes. Both the graph and the image support the localization of the Mouse anti-Human H2Ax Monoclonal - Identifyn™ SRM Alexa Fluor™ 647 IgG1 direct conjugate to the nucleus.

EVIDENCE TIER	Widefield SIM - Apotome
INSTRUMENT	ZEISS Axio Observer.Z1/7 Apotome
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 532; 555; 647
METADATA VALUES	47

STRUCTURED ACQUISITION METADATA / WIDEFIELD SIM - APOTOME

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:
Objective	EC Plan-Neofluar 63X/1.25 oil M27
Z-Stack	86 Slices (8.5 μ m)
Channels	3
Scaling (Per Pixel)	0.143 μ m X 0.143 μ m X 0.100 μ m
Image Size (Pixels)	786 X 689
Image Size (Scaled)	112.29 μ m X 98.43 μ m
Bit Depth	14 Bit
Image Center Position	X: 45.18 mm, Y: 49.96 mm
Stage Position	X: 45.18 mm, Y: 49.96 mm
ROI Center Offset	X: 1.14 μ m, Y: 0.00 μ m
Reflector	50 Cy 5 AF532-49014 49 DAPI
Beam Splitter	660 550 395
Filter Excitation Wavelength	625 - 655 515 - 545 335 - 383
Filter Emission Wavelength	665 - 715 555 - 595 420 - 470
Contrast Method	Fluorescence
Light source	X-Cite Xylis Lamp Intensity @15.00% X-Cite Xylis Lamp Intensity @8.00%
Exposure Time	200 ms 50 ms 20 ms
Excitation Wavelength	653 528 353
Emission Wavelength	668 553 465
Effective NA	1.25
Imaging Device	807 mono
Camera Adapter	1X
Depth of Focus	1.30 μ m 1.07 μ m 0.90 μ m
Binning Mode	2,2
Apotome SIM Image Processing	Display Mod" "Optical Sectioning" - enabled correction selected.
3D Volume Projection	Precise
Appearance	Threshold 13% all channels, Ramp 50% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
Z-Position	47 of 86

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Identifyn® Secondary Antibody

Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat Anti-Rabbit, IgG (H + L)

Cat ID# - SA 000024

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ -H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ -H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ -H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ -H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with Thermo Fisher, Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS.

The second primary antibody, Golgi GOLPH2 Polyclonal Antibody, was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat

Anti-Rabbit, IgG (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI, cat #P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:

Z Stack - (3D)- Panel A

Z-Stack = 86 Slices (8.5 µm)

Channels = 3

Scaling (Per Pixel) = 0.143 µm X 0.143 µm X 0.100 µm

Image Size (Pixels) = 786 X 689

Image Size (Scaled) = 112.29 µm X 98.43 µm

Bit Depth = 14 Bit

Image Center Position = X: 45.18 mm, Y: 49.96 mm

Stage Position = X: 45.18 mm, Y: 49.96 mm

ROI Center Offset = X: 1.14 µm, Y: 0.00 µm

647 -

Reflector = 50 Cy 5

Beam Splitter = 660

Filter Excitation Wavelength = 625 - 655

Filter Emission Wavelength = 665 - 715

Contrast Method = Fluorescence

Light source = X-Cite Xylis Lamp Intensity @15.00%

Exposure Time = 200 ms

Excitation Wavelength = 653

Emission Wavelength = 668

Effective NA = 1.25

Imaging Device = 807 mono

Camera Adapter = 1X

Depth of Focus = 1.30 µm

Binning Mode = 2,2

532 -

Reflector = AF532-49014
 Beam Splitter = 550
 Filter Excitation Wavelength = 515 - 545
 Filter Emission Wavelength = 555 - 595
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @8.00%
 Exposure Time = 50 ms
 Excitation Wavelength = 528
 Emission Wavelength = 553
 Effective NA = 1.25
 Imaging Device = 807 mono
 Camera Adapter = 1X
 Depth of Focus = 1.07 μ m
 Binning Mode = 2,2

DAPI -

Reflector = 49 DAPI
 Beam Splitter = 395
 Filter Excitation Wavelength = 335 - 383
 Filter Emission Wavelength = 420 - 470
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @8.00%
 Exposure Time = 20 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.25
 Imaging Device = 807 mono
 Camera Adapter = 1X
 Depth of Focus = 0.90 μ m
 Binning Mode = 2,2

Post Processing - SIM Processing

Apotome SIM Image Processing = Display Mod" "Optical Sectioning" - enabled correction selected.
 Normalization, Phase Corrector, Fourier Filter, and Deconvolution were NOT used.

3D Image Processing - Panel A

3D Volume Projection = Precise
 Appearance = Threshold 13% all channels, Ramp 50% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Colocalization View - Panel B & C

Z-Position = 47 of 86
 Horizontal Axis - 647, Threshold 824, Range Auto
 Vertical Axis - DAPI, Threshold 279, Range Auto

Summary

Panel A - This panel presents a three-dimensional volume projection of the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647 (red) conjugate and a Golgi primary antibody labeled with Identifyn® SRM Alexa Fluor™ 532 (yellow), utilizing widefield structured illumination microscopy (SIM) resolution.

Panel B - This panel displays a graphical colocalization view featuring the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647 (X-axis) and DAPI (Y-axis).

Panel C - This panel shows the image from which the colocalization measurements were obtained at z-plane 47 out of a total of 86 collected planes. Both the graph and the image support the localization of the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647 IgG1 direct conjugate to the nucleus.

Image Correction

None

Image by E.F. Simon

Image Analysis by B.T. Bennett


Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	Widefield SIM - Apotome
Instrument	ZEISS Axio Observer.Z1/7 Apotome
Microscope configuration	Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:
Structured source metadata values	47
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	3A51AF6482DC5DE5612D357F5D78C3F277B5F0E89924D0520F821D22904EAC55
Method PDF SHA-256	474DB57369216E0162CFB61255FC88938599DDDF72DD6DE07ED0019545720C8F

ORIGINAL IMAGE EVIDENCE / CONFOCAL



EVIDENCE TIER	Confocal
INSTRUMENT	ZEISS LSM 980
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 532; 647
METADATA VALUES	55

STRUCTURED ACQUISITION METADATA / CONFOCAL

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	181 Slices (5.4 μm)
Channels	3
Scaling (Per Pixel)	0.059 μm X 0.059 μm X 0.030 μm
Image Size (Pixels)	1092 X 728
Image Size (Scaled)	64.23 μm X 42.82 μm
Bit Depth	16 Bit
Image Center Position	X: 85.51 mm, Y: 47.62 mm
Stage Position	X: 85.51 mm, Y: 47.62 mm
ROI Center Offset	X: 0.00 μm , Y: 0.00 μm
Reflector	None
Contrast Method	Fluorescence
Pinhole	1.00 AU / 66 μm 1.00 AU / 53 μm 1.00 AU / 42 μm
LSM MBS	MBS 488/639, 405/730
LSM SBS	SBS SP 615 Mirror
Laser Wavelength	639 nm @0.2% 488 nm @0.05% 405 nm @0.2%
Laser Blanking	Enabled
Scan Mode	Frame
Scan Zoom	2.0
Rotation	0°
Sampling	1.20
Pixel Time	0.30 μs
Frame Time	3.69 s
LSM Scan Speed	12
Scan Direction	Bidirectional
Line Step	1
Averaging	8
Excitation Wavelength	653 528 353
Emission Wavelength	668 553 465
Effective NA	1.4
Detection Wavelength	638-745 525-580 410-470
Imaging Device	GaAsP-Pmt3

FIELD	SOURCE-DOCUMENTED VALUE
Detector Type	GaAsP-PMT
Detector Gain	500 V 650 V
Detector Offset	0
Detector Digital Gain	1.0
LSM Plus Mode	7.2 (3D, Auto) 7.8 (3D, Auto) 6.4 (3D, Auto)
3D Maximum Projection	Precise
Appearance	Threshold 1% all channels, Ramp 0% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 100%, Tone Mapping Enabled
Projection	View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Identifyn® Secondary Antibody

Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat Anti-Rabbit, F(ab')₂ (H + L)

Cat ID# - SA 000030

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ-H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ-H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ-H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ-H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with Thermo Fisher, Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS.

The second primary antibody, Golgi GOLPH2 Polyclonal Antibody, was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat

Anti-Rabbit, F(ab')₂ (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:

Single Plane (2D) - Z-plane 79 of 181 total - Panel A

Z Stack - (3D) - Panel A & B

Z-Stack = 181 Slices (5.4 µm)

Channels = 3

Scaling (Per Pixel) = 0.059 µm X 0.059 µm X 0.030 µm

Image Size (Pixels) = 1092 X 728

Image Size (Scaled) = 64.23 µm X 42.82 µm

Bit Depth = 16 Bit

Image Center Position = X: 85.51 mm, Y: 47.62 mm

Stage Position = X: 85.51 mm, Y: 47.62 mm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

647 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 1.00 AU / 66 μ m
 LSM MBS = MBS 488/639, 405/730
 LSM SBS = SBS SP 615
 Laser Wavelength = 639 nm @0.2%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.0
 Rotation = 0°
 Sampling = 1.20
 Pixel Time = 0.30 μ s
 Frame Time = 3.69 s
 LSM Scan Speed = 12
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 8
 Excitation Wavelength = 653
 Emission Wavelength = 668
 Effective NA = 1.4
 Detection Wavelength = 638-745
 Imaging Device = GaAsP-Pmt3
 Detector Type = GaAsP-PMT
 Detector Gain = 500 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 LSM Plus Mode = 7.2 (3D, Auto)

532 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 1.00 AU / 53 μ m
 LSM MBS = MBS 488/639, 405/730
 LSM SBS = Mirror
 Laser Wavelength = 488 nm @0.05%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.0
 Rotation = 0°
 Sampling = 1.20
 Pixel Time = 0.30 μ s
 Frame Time = 3.69 s
 LSM Scan Speed = 12
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 8
 Excitation Wavelength = 528
 Emission Wavelength = 553
 Effective NA = 1.4
 Detection Wavelength = 525-580
 Imaging Device = GaAsP-Pmt3
 Detector Type = GaAsP-PMT
 Detector Gain = 650 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 LSM Plus Mode = 7.8 (3D, Auto)

DAPI -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 1.00 AU / 42 µm
 LSM MBS = MBS 488/639, 405/730
 LSM SBS = Mirror
 Laser Wavelength = 405 nm @0.2%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.0
 Rotation = 0°
 Sampling = 1.20
 Pixel Time = 0.30 µs
 Frame Time = 3.69 s
 LSM Scan Speed = 12
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 8
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Detection Wavelength = 410-470
 Imaging Device = GaAsP-Pmt3
 Detector Type = GaAsP-PMT
 Detector Gain = 500 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 LSM Plus Mode = 6.4 (3D, Auto)

3D Image Processing - Panel B

3D Maximum Projection = Precise
 Appearance = Threshold 1% all channels, Ramp 0% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 100%, Tone Mapping Enabled
 Projection = View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Panel A - 2-dimensional view at z-plane 79 of the 181 captured during image acquisition illustrates the signal and localization of the Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 direct conjugate (red) Also shown as a control for resolution is a known target, Golgi primary antibody labeled with Identifyn® SRM Alexa Fluor™ 532 (yellow).

Panel B - This panel presents a three-dimensional maximum projection of the Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 direct conjugate (red) and the Golgi primary antibody with Identifyn® SRM Alexa Fluor™ 532 (yellow) at confocal resolution.

Panel C - This panel features a clipped image from the Zeiss Zen image acquisition settings for a z-stack comprising

181 planes. The 639 nm laser power is set at 0.2% to achieve optimal imaging of the direct conjugate over a z-stack distance of 3.6 microns. We also note that the Identifyn® SRM Alexa Fluor™ 532 required 0.05%, indicating high brightness for all Identifyn® products via our patent-pending conjugation chemistry.

Both images demonstrate the compartmental localization of the Mouse anti-Human H2AX Monoclonal -Identifyn® SRM Alexa Fluor™ 647, IgG1 direct conjugate within the nucleus, and the punctate visualization we achieve in diffraction-limited confocal resolution. Note that no image edits were made. Transparency in our channel settings shows that very minimal laser power was required, as the direct conjugate was specifically localized and brilliant.

None

Image by J.K. Bennett


Image Analysis by B.T. Bennett

Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE	
Catalog	DC-000035
Stage	Confocal
Instrument	ZEISS LSM 980
Microscope configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Structured source metadata values	55
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	4A2B6D63DCA250A1525B3FC138AE9FB17186D35E182ED3E13F53F72914CDBF07
Method PDF SHA-256	0F6A4C0FDBB685C3CEC6217731150E55C372BD916618359F9A700DB629556427

ORIGINAL IMAGE EVIDENCE / AIRYSCAN



EVIDENCE TIER	Airyscan
INSTRUMENT	ZEISS LSM 980
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 532; 647; 750
METADATA VALUES	63

STRUCTURED ACQUISITION METADATA / AIRYSCAN

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	181 Slices (5.4 μm)
Channels	3
Scaling (Per Pixel)	0.03529 μm X 0.03529 μm X 0.03000 μm
Image Size (Pixels)	1387 X 1105 679 X 653 263 X 278
Image Size (Scaled)	48.94 μm X 38.99 μm 23.96 μm X 23.04 μm 9.28 μm X 9.81 μm
Bit Depth	16 Bit
Image Center Position	X: 84.88 mm, Y: 47.47 mm
Stage Position	X: 84.88 mm, Y: 47.47 mm
ROI Center Offset	X: 0.00 μm , Y: 0.00 μm
Reflector	None
Contrast Method	Fluorescence
Pinhole	5.00 AU / 328 μm 5.00 AU / 265 μm 5.46 AU / 235 μm
LSM MBS	MBS 488/639, 405/730
LSM SBS	SBS LP 640 SBS LP 525 SBS SP 505
Laser Wavelength	639 nm @0.4% 488 nm @0.5% 405 nm @0.2%
Laser Blanking	Enabled
Scan Mode	Frame
Scan Zoom	2.2
Rotation	0°
Sampling	2.00
Pixel Time	1.21 μs
Frame Time	17.06 s
LSM Scan Speed	6
Scan Direction	Bidirectional
Line Step	1
Averaging	4
Excitation Wavelength	653 528 353
Emission Wavelength	668 553 465
Effective NA	1.4
Detection Wavelength	643-735 530-582 350-499
Imaging Device	Airyscan

FIELD	SOURCE-DOCUMENTED VALUE
Detector Type	GaAsP-PMT
Detector Gain	800 V 750 V
Detector Offset	0
Detector Digital Gain	1.0
Airyscan Mode	SuperResolution 9.6 (3D, Auto) SuperResolution 10.0 (3D, Auto) SuperResolution 8.1 (3D, Auto)
3D Surface Projection	Precise
Appearance	Threshold 20%, Ambient Light 0%, Specular Light 100%, Shininess 30% Threshold 2%, Ramp 98%, Maximum 100%
Background	Black
Light	Azimuth 80°, Elongation -130°, Enabled Tone Mapping Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
3D Volume Projection	Precise

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Identifyn® Secondary Antibody

Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat Anti-Rabbit, F(ab')₂ (H + L)

Cat ID# - SA 000030

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ-H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ-H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ-H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ-H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with Thermo Fisher, Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS.

The second primary antibody, Golgi GOLPH2 Polyclonal Antibody, was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat

Anti-Rabbit, F(ab')₂ (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:

Single Plane (2D) - Z-plane 49 of 181 total - Panel A

Z-Stack = 181 Slices (5.4 µm)

Channels = 3

Scaling (Per Pixel) = 35.29 nm X 35.29 nm X 30.00 nm

Image Size (Pixels) = 1387 X 1105

Image Size (Scaled) = 48.94 µm X 38.99 µm

Bit Depth = 16 Bit

Image Center Position = X: 84.88 mm, Y: 47.47 mm

Stage Position = X: 84.88 mm, Y: 47.47 mm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D)- Panel B

Z-Stack = 181 Slices (5.4 µm)

Channels = 3

Scaling (Per Pixel) = 35.29 nm X 35.29 nm X 30.00 nm

Image Size (Pixels) = 679 X 653

Image Size (Scaled) = 23.96 µm X 23.04 µm

Bit Depth = 16 Bit

Image Center Position = X: 84.88 mm, Y: 47.47 mm

Stage Position = X: 84.88 mm, Y: 47.47 mm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D)- Panel D

Z-Stack = 181 Slices (5.4 µm)

Channels = 3

Scaling (Per Pixel) = 35.29 nm X 35.29 nm X 30.00 nm

Image Size (Pixels) = 263 X 278

Image Size (Scaled) = 9.28 µm X 9.81 µm

Bit Depth = 16 Bit

Image Center Position = X: 84.88 mm, Y: 47.47 mm

Stage Position = X: 84.88 mm, Y: 47.47 mm

ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

647 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 5.00 AU / 328 μ m
 LSM MBS = MBS 488/639, 405/730
 LSM SBS = SBS LP 640
 Laser Wavelength = 639 nm @0.4%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.2
 Rotation = 0°
 Sampling = 2.00
 Pixel Time = 1.21 μ s
 Frame Time = 17.06 s
 LSM Scan Speed = 6
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 4
 Excitation Wavelength = 653
 Emission Wavelength = 668
 Effective NA = 1.4
 Detection Wavelength = 643-735
 Imaging Device = Airyscan
 Detector Type = GaAsP-PMT
 Detector Gain = 800 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 Airyscan Mode = SuperResolution 9.6 (3D, Auto)

532 -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 5.00 AU / 265 μ m
 LSM MBS = MBS 488/639, 405/730
 LSM SBS = SBS LP 525
 Laser Wavelength = 488 nm @0.5%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.2
 Rotation = 0°
 Sampling = 2.00
 Pixel Time = 1.21 μ s
 Frame Time = 17.06 s
 LSM Scan Speed = 6
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 4
 Excitation Wavelength = 528
 Emission Wavelength = 553
 Effective NA = 1.4
 Detection Wavelength = 530-582
 Imaging Device = Airyscan
 Detector Type = GaAsP-PMT
 Detector Gain = 800 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 Airyscan Mode = SuperResolution 10.0 (3D, Auto)

DAPI -

Reflector = None
 Contrast Method = Fluorescence
 Pinhole = 5.46 AU / 235 µm
 LSM MBS = MBS 488/639, 405/730
 LSM SBS = SBS SP 505
 Laser Wavelength = 405 nm @0.2%
 Laser Blanking = Enabled
 Scan Mode = Frame
 Scan Zoom = 2.2
 Rotation = 0°
 Sampling = 2.00
 Pixel Time = 1.21 µs
 Frame Time = 17.06 s
 LSM Scan Speed = 6
 Scan Direction = Bidirectional
 Line Step = 1
 Averaging = 4
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.4
 Detection Wavelength = 350-499
 Imaging Device = Airyscan
 Detector Type = GaAsP-PMT
 Detector Gain = 750 V
 Detector Offset = 0
 Detector Digital Gain = 1.0
 Airyscan Mode = SuperResolution 8.1 (3D, Auto)

Post Processing - AiryScan 3D Processing

Display Mode" "SR"
 Super Resolution Auto Filter Selected
 Process 3D Current Image
 Created

Split View - Panel A

The top-left panel shows Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), the top-right panel features Golgi visualized with Identifyn® SRM Alexa Fluor™ 532 Affinity Purified Goat Anti-Rabbit, F(ab')₂ (H + L), the bottom-left panel shows DAPI nuclear staining, and the bottom-right panel is the merged image of all channels. This split-view image was taken at Z-Position 49 of 181.

3D Image Processing - Panel B

3D Surface Projection = Precise
 Appearance = Threshold 20%, Ambient Light 0%, Specular Light 100%, Shininess 30%
 Background = Black
 Light = Azimuth 80°, Elongation -130°, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

3D Image Processing - Panel D

3D Volume Projection = Precise

Appearance = Threshold 2%, Ramp 98%, Maximum 100%

Background = Black

Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 50%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Panel A - Split view Z-plane 49 of the 181 collected illustrates the signal and localization of Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM AlexaFluor™ 647, IgG1 direct conjugate. In the upper left, we see the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 direct conjugate (red). The upper right shows the Golgi primary antibody with Identifyn® SRM Alexa Fluor™ 532 (yellow), while the lower left features DAPI nuclear staining (blue). The lower right presents a merged image of the three channels. A Region of Interest (ROI) is shown to be examined further in Panel B.

Panel B - This panel provides a three-dimensional surface projection of the ROI shown in panel A. The Mouse anti-Human H2Ax Monoclonal - Identifyn™ SRM Alexa Fluor™ 647, IgG1 direct conjugate (red) and the Golgi primary antibody with Identifyn® SRM Alexa Fluor™ 532 (yellow), utilizing super-resolution Airyscan A Region of Interest (ROI) is shown to be examined further in Panel D. Note that no image edits were performed beyond the selection of the surface projection within Zeiss Zen.

Panel C - This panel displays a clipped image from the Zeiss Zen image acquisition settings for this z-stack, which comprises 181 planes. It is noteworthy that the 639 nm laser power is set at 0.4% to achieve optimal imaging of the direct conjugate throughout the z-stack acquisition.

Panel D - This panel shows a three-dimensional volume projection of the ROI shown in panel B. The Mouse anti-Human H2Ax Monoclonal - Identifyn® SRMAlexa Fluor™ 647, IgG1 direct conjugate (red) and the Golgi primary antibody with Identifyn® SRM Alexa Fluor™ 532 (yellow), employing super-resolution Airyscan. Note that no image edits were performed beyond the selection of the surface projection within Zeiss Zen.

For the first time, using the "step down" mode, we observe independent signals and structural detail rather than the pan-nuclear (diffuse) staining seen in diffraction-limited applications such as widefield, widefield SIM, and confocal microscopy. While we do not claim that each independent signal corresponds to a single molecule of H2AX, we emphasize that the enhanced resolution enables the spatial separation of the signals and reveals structural details about H2AX distribution and possibly its connectivity within the nucleus.

Image Corrections

None

Image by J.K. Bennett

Image Analysis by B.T. Bennett


Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	Airyscan
Instrument	ZEISS LSM 980
Microscope configuration	Zeiss LSM 980 Confocal with NIR Detection, Axio Observer.Z1/7, using a Plan-Apochromat 63X/1.40 oil DIC M27, was used with the following parameters for each dye:
Structured source metadata values	63
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	CF12853224E5CDA625FCEC6923C86C18E73EB7FE310A69617A249B552A7E5783
Method PDF SHA-256	F3013E542DCABD5AD021032FF55DA4B2116260A0B1F432E6AC2C3E353C60AE61

ORIGINAL IMAGE EVIDENCE / SIM



EVIDENCE TIER	SIM
INSTRUMENT	ZEISS Elyra
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 647
METADATA VALUES	70

STRUCTURED ACQUISITION METADATA / SIM

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	124 Slices (3.69 μm)
Channels	3
Scaling (Per Pixel)	0.01689 μm X 0.01689 μm X 0.03000 μm
Image Size (Pixels)	5056 X 5056 2010 X 2426
Image Size (Scaled)	85.40 μm X 85.40 μm 33.95 μm X 40.98 μm
Bit Depth	16 Bit
Image Center Position	X: 54.50 mm, Y: 40.07 mm
Stage Position	X: 54.50 mm, Y: 40.07 mm
ROI Center Offset	X: 0.00 μm , Y: 0.00 μm
Reflectors	SR HE LBF 405/488/561/642
Beam Splitter	405, 488, 561, 642
Filter Excitation Wavelength	200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method	Fluorescence
Illumination Wavelength	640 488 405
Channel Name	T1-Axiocam 820-SR T2-Axiocam 820-SR T3-Axiocam 820 #2-SR
Excitation Wavelength	640 488 405
Emission Wavelength	729 525 445
Effective NA	1.4
Detection Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800
Imaging Device	Axiocam 820 Axiocam 820 #2-SR
Camera Adapter	1X
Exposure Time	120 ms 100 ms 50 ms
Depth of Focus	1.13 μm 0.81 μm 0.69 μm
Binning Mode	2,2
Laser Wavelength	640 nm @3.00% 488 nm @0.9% 405 nm @4.00%
Frame Time	1.59 s 1.33 s 676.00 ms
Scan Direction	Unidirectional
Grating	36.5 μm grating
Processing Dimension	3D
Result Sampling	4

FIELD	SOURCE-DOCUMENTED VALUE
Process Sampling	N/A
Algorithm	SIM
SIM Settings	Standard
Auto Sharpness	On
Sharpness	9.8099 (647), 10.2874 (488), 8.5736 (DAPI)
Fitting	Fast Fit
Normalization	Auto
3D Maximum Projection	Precise
Appearance	Threshold 1% all channels, Ramp 0% all channels, Maximum 100% all channels Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels
Background	Black
Light	Brightness 100%, Tone Mapping NOT Enabled Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
Clipping Image Process	Clipping planes are introduced to define regions of special interest.
X-Y	Activate was selected, textured opaque was chosen, and a blue outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.
Position of Clip	-35%
Clip X	0°
Clip Y	0°
3D Volume Projection	Precise

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Identifyn® Secondary Antibody

Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Goat Anti-Rabbit, F(ab')₂ (H + L)

Cat ID# - SA 000018

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ-H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ-H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ-H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ-H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with Thermo Fisher, Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS.

The second primary antibody, Actin Recombinant Rabbit Monoclonal Antibody (JF47-01), was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 488

Affinity Purified Goat Anti-Rabbit, F('b') $\square$ (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.

Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:

Z Stack - (3D) - Panel A

Z-Stack = 124 Slices (3.69 μ m)

Channels = 3

Scaling (Per Pixel) = 16.89 nm X 16.89 nm X 30.00 nm

Image Size (Pixels) = 5056 X 5056

Image Size (Scaled) = 85.40 μ m X 85.40 μ m

Bit Depth = 16 Bit

Image Center Position = X: 54.50 mm, Y: 40.07 mm

Stage Position = X: 54.50 mm, Y: 40.07 mm

ROI Center Offset = X: 0.00 μ m, Y: 0.00 μ m

Z Stack - (3D) - Panel B & C

Z-Stack = 124 Slices (3.69 μ m)

Channels = 3

Scaling (Per Pixel) = 16.89 nm X 16.89 nm X 30.00 nm

Image Size (Pixels) = 2010 X 2426

Image Size (Scaled) = 33.95 μ m X 40.98 μ m

Bit Depth = 16 Bit

Image Center Position = X: 54.50 mm, Y: 40.07 mm

Stage Position = X: 54.50 mm, Y: 40.07 mm

ROI Center Offset = X: 0.00 μ m, Y: 0.00 μ m

Single Plane (2D) - Z-plane 59 of 124 total - Panel A

647 -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 640
 Channel Name = T1-Axiocam 820-SR
 Excitation Wavelength = 640
 Emission Wavelength = 729
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820
 Camera Adapter = 1X
 Exposure Time = 120 ms
 Depth of Focus = 1.13 μm
 Binning Mode = 2,2
 Laser Wavelength = 640 nm @3.00%
 Frame Time = 1.59 s
 Scan Direction = Unidirectional
 Grating = 36.5 μm grating

488 -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 488
 Channel Name = T2-Axiocam 820-SR
 Excitation Wavelength = 488
 Emission Wavelength = 525
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820
 Camera Adapter = 1X
 Exposure Time = 100 ms
 Depth of Focus = 0.81 μm
 Binning Mode = 2,2
 Laser Wavelength = 488 nm @0.9%
 Frame Time = 1.33 s
 Scan Direction = Unidirectional
 Grating = 36.5 μm grating

DAPI -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 405
 Channel Name = T3-Axiocam 820 #2-SR
 Excitation Wavelength = 405
 Emission Wavelength = 445
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820 #2-SR
 Camera Adapter = 1X
 Exposure Time = 50 ms
 Depth of Focus = 0.69 μ m
 Binning Mode = 2,2
 Laser Wavelength = 405 nm @4.00%
 Frame Time = 676.00 ms
 Scan Direction = Unidirectional
 Grating = 36.5 μ m grating

Post Processing - SIM Processing

Processing Dimension = 3D
 Result Sampling = 4
 Process Sampling = N/A
 Channels = 3
 Algorithm = SIM
 SIM Settings = Standard
 Auto Sharpness = On
 Sharpness = 9.8099 (647), 10.2874 (488), 8.5736 (DAPI)
 Fitting = Fast Fit
 Normalization = Auto

3D Image Processing - Panel A

3D Maximum Projection = Precise
 Appearance = Threshold 1% all channels, Ramp 0% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 100%, Tone Mapping NOT Enabled
 Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
 Clipping Image Process = Clipping planes are introduced to define regions of special interest.
 X-Y = Activate was selected, textured opaque was chosen, and a blue outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = -35%

Clip X = 0°

Clip Y = 0°

3D Image Processing - Panel B

3D Volume Projection = Precise

Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels

Background = Black

Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera

Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Panel A - This panel displays a three-dimensional maximum projection of the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 direct conjugate (red), combined with an Actin primary antibody labeled with Identifyn® SRM Alexa Fluor™ 488 (green), and DAPI (blue), utilizing super-resolution SIM imaging. An X-Y cut plane is employed at position-35%, with X and Y at zero-degree angles; the plane is outlined in blue. This cut plane was used to show the nuclear distribution of the H2AX/647 direct conjugate. The image includes a Region of Interest (ROI) denoted by a white box.

Panel B - This panel presents a three-dimensional volume projection from the ROI identified in Panel A, again showing the Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647 (red) and the Actin primary antibody with Identifyn® SRM Alexa Fluor™ 488 (green). It is essential to note that the Zeiss Zen default settings for volume projection were used, with no edits made to the original acquisition. Tone mapping was disabled to ensure the image was not artificially color-mapped or manipulated in any way.

Panel C - This is a two-dimensional image from z-plane 59 of the 124 total images taken to further investigate, in a maximum projection view, or unmodeled structure, that we consistently observe when transitioning from super-resolution Airyscan to SIM, decreasing our axial resolution from 100-120 nm to 60-80 nm in the X and Y dimensions.

This transition reveals additional structural details within the nuclear compartment, particularly regarding the localization of H2AX 647 direct conjugate staining. The findings are strongly supported by the use of a known marker for actin, as SIM has consistently demonstrated the ability to resolve its filamentous structure, as evidenced by these results. Given the control for resolution employing the actin, the organization of the H2AX direct conjugate staining raises intriguing questions about the possible hypothesis of chromatin structural details. However, we make no claims regarding the underlying biology.

Image Corrections

None

Image by P. Raut

Image Analysis by B.T. Bennett


Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	SIM
Instrument	ZEISS Elyra
Microscope configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Structured source metadata values	70
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	22FEDDE5F41FEFA29A90114A5299E49BB2661835AB8B99166EC0C74C22275D48
Method PDF SHA-256	6B14E9E0BCC49C0C94A84F94000BB15E5D699A5A1B9BD449B36DD87DEDB3F164

ORIGINAL IMAGE EVIDENCE / SIM²



EVIDENCE TIER	SIM ²
INSTRUMENT	ZEISS Elyra
CELL MODEL	HeLa
DETECTED DYES	405/DAPI; 488; 647
METADATA VALUES	71

STRUCTURED ACQUISITION METADATA / SIM²

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Objective	Plan-Apochromat 63X/1.40 oil DIC M27
Z-Stack	154 Slices (4.59 µm)
Channels	3
Scaling (Per Pixel)	0.01689 µm X 0.01689 µm X 0.03000 µm
Image Size (Pixels)	5056 X 5056 2290 X 1950 708 X 639
Image Size (Scaled)	85.40 µm X 85.40 µm 38.68 µm X 32.94 µm 11.96 µm X 10.79 µm
Bit Depth	16 Bit
Image Center Position	X: 56.38 mm, Y: 40.74 mm
Stage Position	X: 56.38 mm, Y: 40.74 mm
ROI Center Offset	X: 0.00 µm, Y: 0.00 µm
Reflectors	SR HE LBF 405/488/561/642
Beam Splitter	405, 488, 561, 642
Filter Excitation Wavelength	200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method	Fluorescence
Illumination Wavelength	640 488 405
Channel Name	T1-Axiocam 820-SR T2-Axiocam 820-SR T3-Axiocam 820 #2-SR
Excitation Wavelength	640 488 405
Emission Wavelength	729 525 445
Effective NA	1.4
Detection Wavelength	419 - 470, 503 - 546, 576 - 617, 657 - 800
Imaging Device	Axiocam 820 Axiocam 820 #2
Camera Adapter	1X
Exposure Time	120 ms 100 ms 50 ms
Depth of Focus	1.13 µm 0.81 µm 0.69 µm
Binning Mode	2,2
Laser Wavelength	640 nm @3.00% 488 nm @0.9% 405 nm @7.00%
Frame Time	1.59 s 1.33 s 676.00 ms
Scan Direction	Unidirectional
Grating	36.5 µm grating
Processing Dimension	3D
Result Sampling	4

FIELD	SOURCE-DOCUMENTED VALUE
Process Sampling	4
Algorithm	SIM ²
SIM ² Settings	Standard
Input SNR	Low
Regularization	None
Regularization Weight	0.0000
Iterations	3
Filter	No Filter
Fitting	Fast Fit
Normalization	Auto
3D Volume Projection	Precise
Appearance	Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels Threshold 20%, Ambient Light 0%, Specular Light 100%, Shininess 30%
Background	Black
Light	Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping Azimuth 80°, Elongation -130°, Enabled Tone Mapping
Projection	View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected) View Angle 19°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)
3D Surface Projection	Precise

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7)

Cat ID# - DC 000035

Identifyn® Secondary Antibody

Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Goat Anti-Rabbit, F(ab')₂ (H + L)

Cat ID# - SA 000018

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ-H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ-H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ-H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ-H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with Thermo Fisher, Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS.

The second primary antibody, Actin Recombinant Rabbit Monoclonal Antibody (JF47-01), was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Goat Anti-Rabbit, F(ab')₂ (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000.

Cells were washed 5X in PBS and mounted using ProLong™ Diamond Antifade Mountant with DAPI cat#P36962.
Control Data - Primary Unconjugated Antibody

Microscope

Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:

Single Plane (2D) - Z-plane 78 of 154 total - Panel A

Z-Stack = 154 Slices (4.59 µm)

Channels = 3
Scaling (Per Pixel) = 16.89 nm X 16.89 nm X 30.00 nm
Image Size (Pixels) = 5056 X 5056
Image Size (Scaled) = 85.40 µm X 85.40 µm
Bit Depth = 16 Bit
Image Center Position = X: 56.38 mm, Y: 40.74 mm
Stage Position = X: 56.38 mm, Y: 40.74 mm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D) - Panel B

Z-Stack = 154 Slices (4.59 µm)

Channels = 3
Scaling (Per Pixel) = 16.89 nm X 16.89 nm X 30.00 nm
Image Size (Pixels) = 2290 X 1950
Image Size (Scaled) = 38.68 µm X 32.94 µm
Bit Depth = 16 Bit
Image Center Position = X: 56.38 mm, Y: 40.74 mm
Stage Position = X: 56.38 mm, Y: 40.74 mm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D) - Panel C

Z-Stack = 154 Slices (4.59 µm)

Channels = 3
Scaling (Per Pixel) = 16.89 nm X 16.89 nm X 30.00 nm
Image Size (Pixels) = 708 X 639
Image Size (Scaled) = 11.96 µm X 10.79 µm
Bit Depth = 16 Bit
Image Center Position = X: 56.38 mm, Y: 40.74 mm
Stage Position = X: 56.38 mm, Y: 40.74 mm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

647 -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 640
 Channel Name = T1-Axiocam 820-SR
 Excitation Wavelength = 640
 Emission Wavelength = 729
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820
 Camera Adapter = 1X
 Exposure Time = 120 ms
 Depth of Focus = 1.13 μm
 Binning Mode = 2,2
 Laser Wavelength = 640 nm @3.00%
 Frame Time = 1.59 s
 Scan Direction = Unidirectional
 Grating = 36.5 μm grating

488 -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 488
 Channel Name = T2-Axiocam 820-SR
 Excitation Wavelength = 488
 Emission Wavelength = 525
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820
 Camera Adapter = 1X
 Exposure Time = 100 ms
 Depth of Focus = 0.81 μm
 Binning Mode = 2,2
 Laser Wavelength = 488 nm @0.9%
 Frame Time = 1.33 s
 Scan Direction = Unidirectional
 Grating = 36.5 μm grating

DAPI -

Reflectors = SR HE LBF 405/488/561/642
 Beam Splitter = 405, 488, 561, 642
 Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
 Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Contrast Method = Fluorescence
 Illumination Wavelength = 405
 Channel Name = T3-Axiocam 820 #2-SR
 Excitation Wavelength = 405
 Emission Wavelength = 445
 Effective NA = 1.4
 Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
 Imaging Device = Axiocam 820 #2
 Camera Adapter = 1X
 Exposure Time = 50 ms
 Depth of Focus = 0.69 μ m
 Binning Mode = 2,2
 Laser Wavelength = 405 nm @7.00%
 Frame Time = 676.00 ms
 Scan Direction = Unidirectional
 Grating = 36.5 μ m grating

Post Processing - SIM² Processing

Processing Dimension = 3D
 Result Sampling = 4
 Process Sampling = 4
 Channels = 3
 Algorithm = SIM²
 SIM² Settings = Standard
 Input SNR = Low
 Regularization = None
 Regularization Weight = 0.0000
 Iterations = 3
 Filter = No Filter
 Fitting = Fast Fit
 Normalization = Auto

Split View - Panel A

The top-left panel features Mouse anti-Human H2Ax Monoclonal - Identifyn® SRM Alexa Fluor™ 647, IgG1 - (1H1D9/B7), the top-right panel shows Actin visualized with Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Goat Anti-Rabbit, F(a") $\square$ (H + L), the bottom-left panel features DAPI nuclear staining, and the bottom-right panel features the merged image of all channels. This image was taken at Z-Position 78 of 154.

3D Image Processing - Panel B

3D Volume Projection = Precise

Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels

Background = Black

Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera

Separation 5% (deselected), Parallax Shift 0% (deselected)

3D Image Processing - Panel C

3D Surface Projection = Precise

Appearance = Threshold 20%, Ambient Light 0%, Specular Light 100%, Shininess 30%

Background = Black

Light = Azimuth 80°, Elongation -130°, Enabled Tone Mapping

Projection = View Angle 19°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera

Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Panel A - This panel illustrates a split view demonstrating the signal and localization of the Mouse anti-Human H2AX Monoclonal -Identifyn® SRM Alexa Fluor™ 647, an IgG1 direct conjugate. The upper left section displays the H2AX direct conjugate (red), while the upper right shows the Actin primary antibody labeled with Identifyn® SRM Alexa Fluor™ 488 (green). The lower left features DAPI nuclear staining (blue), and the lower right presents a merged image of all three channels. A Region of Interest is indicated by the white box, which will be explored further in Panel B.

Panel B - This panel presents a three-dimensional volume projection of the Region of Interest from Panel A. It includes the Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647 (red) and the Actin primary antibody with Identifyn® SRM Alexa Fluor™ 488 (green), and pseudo colored white for visual contrast. This was achieved using super-resolution SIM Squared, with no signal manipulation beyond selecting the volume view in Zeiss Zen software. A Region of Interest is again indicated by the white box, which will be examined further in Panel C.

Panel C - This panel presents a three-dimensional surface projection of the Region of Interest from Panel B. It includes the Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647, which has been pseudo-colored to light grey in Zeiss Zen software to enhance visualization of the structural details against a black background, utilizing super-resolution SIM Squared.

It is noteworthy that our observations with the H2AX mouse monoclonal primary antibody, counterstained with a secondary fluorescent antibody, as well as the other direct conjugates across all samples and images employing super-resolution techniques such as Airyscan, SIM, SIM Squared, and STORM, show consistent structural details. This consistency validates the three-dimensional localization and structural integrity of the images.

While this image may suggest chromatin organization due to the ubiquitous nature of H2AX in the nucleus, this interpretation is hypothetical. All images demonstrate this pattern, which appears to be biologically relevant, but it is not conclusive. Nonetheless, it could serve as a valuable tool in colocalization studies with known protein targets at this level of resolution, approximately 40 nm in the X and Y dimensions.

Image Corrections

Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Goat Anti-Rabbit, F(a'')₂ (H + L) was pseudo-colored in Zen Software from green, the default color, to light grey to provide visual contrast between the 647/488 channels.

Image by P. Raut

Image Analysis by B.T. Bennett


Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	SIM ²
Instrument	ZEISS Elyra
Microscope configuration	Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:
Structured source metadata values	71
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	E3DC894B9946EA60699A815B2E563520EE97EC2C593B08DD8721DA579F9CE7E0
Method PDF SHA-256	194F54AD0383741F22D034D8A87434B64FC066440B727F58609E6F7689EFA9D9

ORIGINAL IMAGE EVIDENCE / SINGLE-MOLECULE STORM



EVIDENCE TIER	Single-molecule STORM
INSTRUMENT	Bruker Vutara VXL
CELL MODEL	HeLa
DETECTED DYES	647
METADATA VALUES	80

STRUCTURED ACQUISITION METADATA / SINGLE-MOLECULE STORM

The values below were deterministically extracted from the preserved Identifyn method PDF text. They remain source-documented metadata, not new BENNETT measurements.

FIELD	SOURCE-DOCUMENTED VALUE
Platform / configuration	Bruker Vutara VXL Biplane Single Molecule Microscope, using an Olympus 60X/1.30 silicon oil objective, was used with the following parameters:
Objective	Olympus 60X/1.30 silicon oil objective
Number of Probes	1
Number of Simultaneous Probes	1
Number of Reference Probes	1
Camera	Both Cameras
Color Channel	Both Channels
Probe Capture Mode	Interleaved
Z -Stack Capture Mode	Interleaved
Multiple Location Capture	Not selected
Post Record Reference	Not Selected
Fluidics	Not Selected
Autofocus Drift Correction	On
Biplane Module	635nm Dichroic
Objective	Silicon Selected
Stage Insert	Rectangular Selected
Course Focus Position	8-Well Chamber Selected
Dichroic	SuperRes
PSF	(Red, Orange)
Heater	Current: 26
Recording Vertical Field of View	100%
Widefield Camera Exposure Time	100ms 20ms
Super Resolution Camera Exposure Time	20ms
Widefield Camera Binning	1
Per-Probe Exposure Time	Selected
Predefined Probe	First probe: 640 nm (635 nm Dichroic); First reference probe: 640 nm (635 nm Dichroic)
Photoactivation/Imaging Laser	0% Selected
Exposure Time	20ms Selected
Record Workflow	Selected Options
View Mode	Widefield 200µm
Capture Mode	Live
Piezo Current	179.2 µm

FIELD	SOURCE-DOCUMENTED VALUE
Widefield 640nm	0.1X Neutral Density Filter: Selected, @0.2% Laser Input
Reference Probe 1 Laser control (640 nm (640 nm Dichroic))	Not Selected
Widefield camera Exposure Time	100ms
Widefield Laser	"Widefield 640nm" was deselected, and "calibrate" under Autofocus was selected.
Current Intensity	4.5
Calibration Value	-5070 nm/unit
Autofocus Mode	Pre-Capture
SuperRes	50µm is Selected
Piezo Record	Begin: 178; End: 178.6
SuperRes 640nm	0.1X Neutral Density Filter: Selected, @0.02% Laser Input
Record Workflow - Record Tab	Capture Parameters
Frames	Probe 1: 300
Z Positions	3 using Steps: 0.2µm (200nm)
Cycles	20
Time Points	1
Pause Between Cycles	Not Selected
Automatic Photoactivation Laser Control	Not Selected
Laser Power	30%
Cutout Width	16px
Max Frame Accumulation	1 frame
Max Accumulation Offset	1px
Fiducial at Max Accumulation	Not selected
Max Particles Per Frame	No Max
Background SLPE Removal	Not Selected
Cutout/Background Removal	Not Selected
Preview Cutout Positions	Selected
Probe1	BG threshold: 7
Sizing Mode	Radial Precision
Particle size	4nm
Color by	Constant
Color Table	Single
Opacity Mode	Constant; Opacity 1: 03
Front-to-back Attenuation	Not Selected
Method	Particle (direct)
Time Intervals	7 Intervals

FIELD	SOURCE-DOCUMENTED VALUE
Pixel Size	50nm
Smoothing	8.00
PSF Z offset	-600 to +600
Radial precision	30
Axial precision	70
Clustering Algorithm	DBSCAN
Maximum Particle Distance	0.1µm (1000nm)
Maximum Particle Count to Form Cluster	10
Computer Hulls	Selected
Compute Density Isosurfaces	Not Selected
Edit Settings	All Probes
Hull Alpha shape Radius	1000nm
Compute	Selected
Following Cluster Analysis Computation, under the Results Tab, ID Unassigned	Deselected

SOURCE METHOD

The text below is carried from the preserved Identifyn source record. Typography and pagination have been normalized for this preservation dossier; scientific wording has not been editorially rewritten.

Identifyn Standards for Optical Microscopy Methods™

Identifyn Primary Antibody

Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647 - Direct Conjugate

Cat ID# - DC 000035

H2AX, a variant of the histone protein H2A, and its primary function relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ -H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ -H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ -H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ -H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown in a μ -Slide 8 Well high Glass-Bottom Ibidi chamber. DNA damage was induced via treatment with 2 μ M Etoposide-spiked complete growth media for 6 hours, followed by -20°C methanol fixation and blocked with Thermo Fisher Fish Serum Blocker cat#37527X3. The primary antibody, Mouse anti-Human H2AX Monoclonal - Identifyn® SRM Alexa Fluor™ 647 - Direct Conjugate, was incubated for 1 hour at room temperature at a dilution of 1:100. Cells were washed 5X in Identifyn's™ proprietary Wash Buffer, followed by the addition of Identifyn's™ proprietary STORM Buffer to the sample chamber, enabling image acquisition.

Microscope

Bruker Vutara VXL Biplane Single Molecule Microscope, using an Olympus 60X/1.30 silicon oil objective, was used with the following parameters:

Calibration -

Calibration follows Identifyn's Calibration Standard Operating Procedure, which can be found here.

<https://identifyn.com/storm-calibration>

Image Acquisition -

SRX Software - Under the 'Single Molecule Localization' tab, 'Record' is selected and contains the following workflow tabs and parameters.

Configuration Tab - Capture Configuration

Number of Probes: 1
 Number of Simultaneous Probes: 1
 Number of Reference Probes: 1
 Camera: Both Cameras
 Color Channel: Both Channels
 Probe Capture Mode: Interleaved
 Z -Stack Capture Mode: Interleaved
 Multiple Location Capture: Not selected
 Post Record Reference: Not Selected
 Fluidics: Not Selected
 Autofocus Drift Correction: On

Microscope configuration

Biplane Module: 635nm Dichroic
 Objective: Silicon Selected
 Stage Insert: Rectangular Selected
 Course Focus Position: 8-Well Chamber Selected
 Dichroic: SuperRes
 PSF: (Red, Orange)
 Heater: Current: 26

Camera Configuration

Recording Vertical Field of View: 100%
 Widefield Camera Exposure Time: 100ms
 Super Resolution Camera Exposure Time: 20ms
 Widefield Camera Binning: 1
 Per-Probe Exposure Time: Selected

Probe Configuration

Predefined Probe: First probe: 640 nm (635 nm Dichroic); First reference probe: 640 nm (635 nm Dichroic)
 Photoactivation/Imaging Laser: 0% Selected
 Exposure Time: 20ms Selected

Experiment Configuration

Record Workflow: Selected Options
 View Mode: Widefield 200µm
 Capture Mode: Live

Stage and Piezo Selection

Examination of Sample

Cell(s) or Cellular Structure Localization - The red-bordered box identifies the area of the SuperRes view and defines a suitable cell(s) or cellular structure. The cell(s) or structures are optimally focused.

Piezo Current: 179.2 μ m

Probe 1 Laser control (640 nm (635 nm Dichroic)

Widefield 640nm: 0.1X Neutral Density Filter: Selected, @0.2% Laser Input
Reference Probe 1 Laser control (640 nm (640 nm Dichroic): Not Selected

Advanced Tools Option

Widefield camera Exposure Time: 100ms
Default values are maintained unless otherwise noted.

Widefield Laser: "Widefield 640nm" was deselected, and "calibrate" under Autofocus was selected.

Autofocus - Calibration Values

Current Position" -0.979
Current Intensity: 4.5
Calibration Value: -5070 nm/unit
Autofocus Mode: Pre-Capture

View Mode

SuperRes: 50 μ m is Selected

Capture Mode

Live

Stage and Piezo Selection

Cell(s) or Cellular Structure Localization - Positioning of the target Cell or Subsection of the Target Cell is optimized and optimally focused.

Piezo Record: Begin: 178; End: 178.6

Probe 1 Laser control (640 nm (635 nm Dichroic)

SuperRes 640nm: 0.1X Neutral Density Filter: Selected, @0.02% Laser Input

Reference Probe 1 Laser control (640 nm (635 nm Dichroic)

Not Selected

Advanced Tools

Widefield Camera Exposure Time: 20ms
Default values are maintained unless otherwise noted.

Record Workflow - Record Tab: Capture Parameters

Frames: Probe 1: 300

Z Positions: 3 using Steps: 0.2µm (200nm)

Cycles: 20

Time Points: 1

Pause Between Cycles: Not Selected

Automatic Photoactivation Laser Control: Not Selected

Probe 1 Laser control (640 nm (640 nm Dichroic)

Laser Power: 30%

Expert Configuration Options - The Localization Tab

Beyond the values shown, all other parameters were set to their default values.

Localization

Cutout Width: 16px

Max Frame Accumulation: 1 frame

Max Accumulation Offset: 1px

Fiducial at Max Accumulation: Not selected

Max Particles Per Frame: No Max

Background SLPE Removal: Not Selected

Cutout/Background Removal: Not Selected

Preview Cutout Positions: Selected

Record Workflow - Localization Tab

Data Collection

The folder is created automatically with the localization file for this image acquisition.

Recording Summary

All parameters were carried over from previous tabs and are kept except:

Probe1: BG threshold: 7

Region of Interest -

Not selected

Display

Was deployed to inspect if the BG threshold applied chooses localizations in an expected manner

Record Workflow - View Tab

Region of Interest

Default values

Visualization Parameters

Point Splatting

Sizing Mode: Radial Precision

Particle size: 4nm

Color by: Constant

Color Table: Single

Opacity Mode: Constant; Opacity 1: 03

Front-to-back Attenuation: Not Selected

Post Acquisition analysis - View Tab

Drift Correction - Computed with the following parameters:

Method: Particle (direct)

Time Intervals: 7 Intervals

Pixel Size: 50nm

Smoothing: 8.00

Advanced Particle Filters - Applied to the drift-corrected data with the following parameters

PSF Z offset: -600 to +600

Radial precision: 30

Axial precision: 70

The other parameters were not applied.

Statistical analysis - Computed to clean low background signals.

Statistics Dialog - Cluster Analysis

Clustering Algorithm: DBSCAN

Parameters

Maximum Particle Distance: 0.1 μ m (1000nm)

Maximum Particle Count to Form Cluster: 10

Computer Hulls: Selected

Compute Density Isosurfaces: Not Selected

Edit Settings: All Probes

Hull Alpha shape Radius: 1000nm

Compute: Selected

Following Cluster Analysis Computation, under the Results Tab, ID Unassigned: Deselected

Image by S. Thakur

© Copyright 2026 GMD12 LLC.

SOURCE PROVENANCE

Catalog	DC-000035
Stage	Single-molecule STORM
Instrument	Bruker Vutara VXL
Microscope configuration	Bruker Vutara VXL Biplane Single Molecule Microscope, using an Olympus 60X/1.30 silicon oil objective, was used with the following parameters:
Structured source metadata values	80
Metadata source	PRESERVED_SOURCE_METHOD_PDF_TEXT
Cell model	HeLa
Image SHA-256	848AEBF576CE9A02764A75EC2B616486D08E8D607D517B42C43704B2208109D6
Method PDF SHA-256	D82DDE40727F30E4C81CEA5FEC008924D5AFD4A022CB2BE9ECD4948169DBECF2